Transcutaneous Electrical Nerve Stimulation (TENS)

Model GF-3 / GF-3T

Operation Manual

Read this manual before operating your GF-3 / GF-3T.

Save this manual for future use.

The most current version of this manual can be found online at
www.grahamfield.com

CONTENTS

GENERAL DESCRIPTION	3
WHAT IS TENS?	3
INDICATIONS AND CONTRAINDICATIONS	3
SAFETY	4
ABOUT THE DEVICE.....	6
EXPLANATION OF CONTROL FUNCTIONS	7
ATTACHING THE LEAD WIRES	8
ELECTRODE SELECTION AND CARE.....	8
TIPS FOR SKIN CARE	8
CONNECTING THE DEVICE	9
BATTERY INFORMATION	10
CARING FOR YOUR DEVICE.....	11
TROUBLESHOOTING	11
TECHNICAL SPECIFICATIONS.....	12
OUTPUT SPECIFICATIONS.....	12
LIMITED WARRANTY	13
MANUEL D'UTILISATION EN FRANÇAIS.....	17

GF Health Products, Inc. is not responsible for typographical errors. Packaging, warranties, products and specifications are subject to change without notice.

Graham-Field and **Grafcō** are registered trademarks of GF Health Products, Inc.

GENERAL DESCRIPTION

TENS, Transcutaneous Electrical Nerve Stimulation, is a method of relieving symptomatic chronic intractable pain.

WHAT IS TENS?

TENS is a treatment whereby electrical impulses are applied to nerves through electrode pads placed on the skin. TENS is non-invasive and does not use pharmaceuticals.

TENS uses a two-pronged approach to pain relief. First, sensory nerves are targeted, stimulating them to block pain signals and prevent their transmission to the brain. Second, TENS promotes the production of endorphins — neurochemicals occurring naturally in the brain — which have analgesic properties.

INDICATIONS AND CONTRAINDICATIONS

Read the operation manual before using this TENS device.

Federal law (USA) restricts this device to sale by or on the order of a physician.

Indications

TENS is indicated to be used under a physician's prescription for the symptomatic relief of chronic intractable pain.

Contraindications

- Any electrode placement which applies current to the carotid (neck) region.
- Patients with implanted electronic devices (for example, a pacemaker) or metallic implants should not undergo TENS treatment without first consulting a physician.
- Any electrode placement which causes current to flow transcerebrally (through the head).
- The use of unit whenever pain symptoms are undiagnosed and the etiology is unknown.

SAFETY

Always follow basic safety precautions, including the following:

- ⚠ WARNING: Indicates a potential hazard situation or unsafe practice that, if not avoided, could result in death or serious personal injury.**
- ▲ CAUTION: Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury or product / property damage.**

Warnings

- ⚠ WARNING: Explosion hazard**
Explosion hazard is possible if used in the presence of explosives, flammable materials or flammable anesthetics.
- ⚠ WARNING: Heart disease**
Caution should be used when applying the device to patients suspected of having heart disease. Further clinical data is needed to show if there are adverse side effects on individuals with heart disease.
- ⚠ WARNING: Keep this device out of the reach of children.**
- ⚠ WARNING: The safety of the device during pregnancy or delivery has not been established.**
- ⚠ WARNING: Do not place electrodes on front of the throat. This may result in spasms of the laryngeal and pharyngeal muscles.**
- ⚠ WARNING: Do not place the electrodes over the carotid nerve.**
- ⚠ WARNING: The device is not effective for pain of central origin (headaches).**
- ⚠ WARNING: Avoid adjusting controls while operating machinery or vehicles.**
- ⚠ WARNING: The device may interfere with electronic monitoring equipment (such as ECG monitors and ECG alarms).**
- ⚠ WARNING: Do not change any mode during treatment.**

- ⚠ WARNING: Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. Using an alternate conductive medium or alternate electrode placement can usually reduce the irritation. Consult your physician / clinician before using an alternative conductive medium or electrode placement.**
- ⚠ WARNING: Electrodes should not be placed over the eyes, in the mouth, or internally.**
- ⚠ WARNING: The device has no curative value.**
- ⚠ WARNING: TENS devices should be used only under the continued supervision of a physician / clinician.**
- ⚠ WARNING: TENS is a symptomatic treatment and as such suppresses the sensation of pain which would otherwise serve as a protective mechanism.**

Precautions / Adverse Reactions

- ▲ CAUTION: Isolated cases of skin irritation may occur at the site of electrode placement following long-term application.**
- ▲ CAUTION: If skin irritation occurs TENS treatment should be stopped and electrodes removed until the cause of the irritation can be determined.**
- ▲ CAUTION: Effectiveness is highly dependent upon patient selection of a doctor qualified in the management of pain patients.**
- ▲ CAUTION: If the device treatment becomes ineffective or unpleasant, stimulation should be discontinued until reevaluation by a physician / clinician.**
- ▲ CAUTION: Always turn the device OFF before applying or removing electrodes.**

ABOUT THE DEVICE


This device is a battery-operated device that includes two controllable output channels. This device creates electrical impulses in which intensity, duration, and modulation can be altered. The device controls are easy to use and the slide cover protects accidental changes in settings.

System Components

Your device will include the following components or accessories:

- TENS Unit
- Carrying case
- Lead wires
- 9-Volt battery
- Operation Manual
- Electrodes

Device Controls

Front Panel Cover

This cover located on the front of the unit conceals the controls for Pulse Width, Pulse Rate, Mode Selector, and Timer (GF-3T only). Press the top portion of the cover and pull down in order to open the cover.

EXPLANATION OF CONTROL FUNCTIONS

Ch1 / Ch2 Knobs	The Channel 1 and Channel 2 Intensity Control Knobs control the strength of the stimulation and also function as ON/OFF controls for each channel.	
MODE Selection Switch	The MODE Selection Switch sets the treatment mode. The three modes are <u>Burst (B)</u> , <u>Continuous (C)</u> , and <u>Modulation (M)</u> .	<u>Burst Mode (B)</u> releases individual bursts twice per second, pulse width is adjustable, and the pulse rate is set at 100 Hz per second.
		<u>Continuous Mode (C)</u> delivers stimulation continuously at the settings determined by the intensity, rate, and width knobs.
		<u>Modulation Mode (M)</u> decreases the pulse width to 60% of the original setting; this decreased pulse width is maintained for 1.5 seconds before returning to the original pulse width setting, which is maintained for 3.5 seconds. The cycle is then repeated. The intensity and pulse rate are adjustable.
Pulse Width Control Knob	The Pulse Width Control Knob regulates the pulse width for both channels.	
Pulse Rate Control Knob	The Pulse Rate Control Knob regulates the number of pulses per second for both channels.	
Timer Selection Switch * (GF-3T Only)	The Timer Selection Switch selects the duration of TENS treatment. The three Timer selections are <u>Continuous (C)</u> , <u>30 minutes (30)</u> , and <u>60 minutes (60)</u> . * To resume operation or reset the Timer, turn the Intensity Control Knob OFF and then ON.	

ATTACHING THE LEAD WIRES

⚠ WARNING: Ensure the device is OFF before connecting the lead wires.

⚠ WARNING: Never insert the lead wire plug into an AC power supply socket. Personal injury and/or damage to the TENS unit could occur.

▲ CAUTION: Use care when you plug and unplug the wires. Pulling on the lead wire instead of its insulated connector may cause wire breakage.

The lead wires provided with the device insert into the ports located the top of the unit. After connecting the wires to the unit, attach each wire to an electrode.

Lead wires provided with the device are compliant with mandatory compliance standards set forth by the FDA.

ELECTRODE SELECTION AND CARE

Using Electrodes

Use the electrodes as prescribed. Follow application procedures outlined in electrode packaging to maintain stimulation and prevent skin irritation. The electrode packaging provides instructions for care, maintenance, and proper storage of electrodes.

TIPS FOR SKIN CARE

Good skin preparation is important for effective and comfortable use of your TENS device.

- Always clean the electrode site with mild soap and water solution, rinse well, and dry thoroughly prior to any electrode application.
- Any excess hair should be clipped, not shaved, to ensure good electrode contact with the skin.
- If a skin treatment or preparation is recommended by your physician / clinician, apply the skin treatment as recommended, let dry, and apply electrodes as directed. Following these recommendations will both reduce the chance of skin irritation and extend the life of your electrodes.

- Avoid excessive stretching of the skin when applying electrodes. Proper application is best accomplished by applying the electrode, then smoothly pressing it in place from the center outward.
- When removing electrodes, always remove by pulling in the direction of hair growth.
- It may be helpful to rub skin lotion on electrode placement area when not wearing electrodes.

CONNECTING THE DEVICE

Insert battery

Turn the device to the OFF position before inserting or removing the battery. When inserting the battery, ensure the battery polarity (+ and -) markings match the markings on the device.

Prepare the Skin

Prepare the skin as previously described and according to the instructions provided with your electrodes. Before attaching the electrodes, identify the area that your physician / clinician has recommended for electrode placement.

1. Connect the Lead Wires to the electrodes: Connect the lead wires to the electrodes before applying the electrodes to the skin.

⚠ WARNING: Ensure both intensity controls for Channel 1 and 2 are turned to the “OFF” Position (counterclockwise) before applying the electrodes.

2. Place Electrodes on the Skin: Place the electrodes on the skin as recommended by your physician / clinician.
3. Insert Lead Wire Connector into the Device: Plug end of lead wire into the channel output port (jack) to be used; push the plug in as far as it will go.
4. Select Treatment Settings: Ensure your unit is still set to the proper settings recommended by your physician / clinician.

5. **Adjusting Channel Intensity Control:** Locate the intensity control knob (Channel 1 or 2) at the top of the unit. Slowly turn the intensity control knob clockwise until the stimulation is at the level recommended by your physician / clinician (if you don't feel anything, turn the knob OFF then ON again and carefully turn the control knob until you feel a tingling or slight twitch under or around the electrodes). Always start with the lowest setting and increase the intensity slowly.

If the stimulation levels are uncomfortable or become uncomfortable, reduce the stimulation intensity to a comfortable level; or cease stimulation and contact your physician.

BATTERY INFORMATION

A 9-volt battery is provided with your unit. When the indicator lights on the front of the unit no longer illuminate, the battery has become too weak to power the unit, and the existing battery should be replaced with a new battery. At this point, the unit will turn OFF until a new battery is inserted.

Replacing the Battery

When the indicator lights on the front of the unit no longer illuminate, the battery should be replaced.

1. Turn unit OFF.
2. Remove the front panel cover by pressing on the top of the panel and pressing down in order to slide the panel down. Continue sliding the panel downwards until the panel is completely removed from the unit. This will reveal the battery compartment.
3. Remove the discharged battery from the device.
4. Place new battery in the compartment. Note: Be sure the proper polarity (+ and -) markings match the markings in the device.
5. Replace the front panel cover.
6. Dispose of the old battery according to local guidelines and regulations.

CARING FOR YOUR DEVICE

Your device may be cleaned by wiping gently with a damp cloth moistened with mild soap and water. Do not immerse the device in water or other liquids.

Wipe lead wires with damp cloth moistened with soap and water. Do not immerse the lead wires.

To properly store the device for an extended period of time, remove the battery from the unit. Place the unit and accessories in the carrying case provided and store in a cool, dry location.

TROUBLESHOOTING

If the device does not function properly:

1. Ensure the battery is properly installed or replace the battery. Be sure to observe proper polarity markings when replacing the battery. If the indicator lights on the front of the unit no longer illuminate when it is turned on, replace the battery and check again.
2. If the indicator light is flashing, the intensity has been adjusted, and no stimulation is felt, check to ensure the lead wires are properly connected and the electrodes are properly applied to the skin. If the unit appears to be functioning and no stimulation is felt, the lead wires or electrodes may need to be replaced.
3. If the battery appears to be charged and the unit is not functioning, turn both intensity control knobs to the OFF position (counterclockwise). Then gradually turn the intensity Control Knob (clockwise) until stimulation is felt. If device still is not working, turn the unit off and contact your authorized GF Health Products, Inc distributor.

If there is any other problem, please contact an authorized GF Health Product, Inc. distributor. Do not try to repair a defective device.

TECHNICAL SPECIFICATIONS

Channel	Dual, isolated between channels
Modes of Operation	Burst, Continuous, and Modulation
Pulse Intensity	Adjustable 0 mA – 80 mA peak into 500 ohm load each channel, constant current
Pulse Rate	2 Hz-150 Hz (adjustable)
Pulse Width	30uS-260uS (adjustable)
Timer (GF-3T only)	Continuous, 30 minutes, and 60 minutes
Burst Mode (B)	Burst consists of 2 bursts per second at 100 Hz
Wave Form	Asymmetrical Bi-Phasic Square Pulse
Voltage	0 - 100 Volt (open circuit)
Power Source	9-Volt battery
Dimensions	95 mm x 65 mm x 23.5 (Height x Width x Thickness)
Poids	115 g (including battery)

OUTPUT SPECIFICATIONS

Mode	Intensity (mA)	Pulse Width (uSec)	Pulse Rate Frequency (Hz)	Cycle Time (Sec)
Burst	Adj. 0-80	Adj. 30-260	100 Hz fixed N/A 2 bursts per Sec	
Continuous	Adj. 0-80	Adj. 30-260	Adj. 2-150 Hz	N/A
Modulation	Adj. 0-80	Modulates down from preset width setting by 60% then back to original setting	Adj. 2-150 Hz	5 Sec total time

LIMITED WARRANTY

GF Health Products, Inc. warrants the Transcutaneous Electrical Nerve Stimulation (TENS), Model GF-3 / GF-3T, against manufacturer's defects for one year.

The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alterations or disassembly by unauthorized individuals.

During the warranty period, defective items will be repaired or replaced at GF Health Products, Inc. option. Warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to GF Health Products, Inc.

USA Corporate Headquarters:

GF Health Products, Inc.
2935 Northeast Parkway
Atlanta, Georgia 30360 USA

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© 2010 GF Health Products, Inc.