

10-year Calibration Warranty

This blood pressure unit is guaranteed to remain accurate for 10 years. This warranty entitles the original purchaser to readjustment, if the unit becomes inaccurate under normal conditions of use. Readjustment does not include the replacement of any broken parts or any part of the inflation system.

GF Health Products, Inc.
2935 Northeast Parkway
Atlanta, GA 30360
1-800-347-5678



A Graham-Field Brand

How To Use Your

Latex Free

Labtron

Sphygmomanometer

GF1000096RevA10

Manometer Gauge:

Your sphygmomanometer is equipped with a 300mmHg aneroid manometer. Requiring minimal maintenance, the gauge is guaranteed accurate to ± 3 mmHg over its full range. The gauge may be cleaned with a slightly damp cloth, but should not be opened or dismantled under any circumstances. Gauge accuracy can be checked visually; simply be certain the needle rests in the printed oval when the unit is fully deflated.

10-Year Calibration Warranty:

Units subject to negligence, misuse, improper operation, improper maintenance or improper storage, normal wear, failure to follow instructions or damages beyond GF Health Products, Inc.'s control are not covered by the warranty, and that evaluation shall be determined solely by GF Health Products, Inc.

In the event that the unit is in need of calibration, carefully wrap the unit, include a note containing your name, return address, daytime phone number, a description of the malfunction, and proof of purchase (copy of receipt), and return it to:

GF Health Products, Inc.
2935 Northeast Parkway
Atlanta, GA 30360

For a shipping and handling charge of \$7.00, the gauge will be recalibrated. If you have questions, please call Graham-Field Customer Service at 1-800-347-5678.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

This unit is the result of extensive research, and uses only the finest, most durable materials available. It is designed to provide you with many years of accurate, dependable service with minimal care and maintenance.

Your Labtron sphygmomanometer consists of an aneroid manometer (gauge); complete inflation system (includes cuff, inflation bladder, squeeze bulb, and air release valve); zippered carrying case; and a 10-year calibration warranty.

Applying the Cuff:

- Place the cuff over bare upper arm, just above the bend of the inner elbow.
- Wrap the cuff snugly around arm and engage hook & loop fastener.
- To verify a correct fit, ensure that two fingers fit between the cuff and skin surface. If sufficient room is not available, readjust the cuff. A different size may be needed if readjustment does not provide sufficient room or if the cuff cannot be secured tightly enough.
- Please note that these cuffs are pregauged for proper arm size. If the index line falls outside the specified range, please ensure you are using the proper size cuff (i.e. adult, infant, etc.).

Overinflation:

Inflating the instrument beyond 300mmHg can result in damage to manometer and inflation system.

Storage:

Wrap cuff around manometer and bulb; store in zippered carrying case provided.

Care & Maintenance:

The cuff can be easily hand washed. However, remove the bladder and ensure hook & loop are fastened before washing. Do not iron fasteners.