

# How To Replace Your Lumex® Set N' Go Rollator Brakes



1. The following items are needed to replace your handbrakes:

- Channel Lock Pliers
- Thick Towel or Rag
- Replacement Handbrake

2. Wrap towel around handgrip for protection, then use pliers to apply pressure and rotate handgrip down as shown. You will hear a pop as handgrip adhesive releases. Continue twisting to slide handgrip off.

You can now remove the top screw shown above, that secures the brake assembly onto the handle tubing.


3. You will see a spring and cable mechanism on the rear wheels. Identify the brass cable lock and loosen the socket head cap screw several turns; you can then remove the brass cable lock by pulling straight downward. Pull firmly to remove the end cable cap. You will have the pieces shown in the right side picture.

4. Install the brake handle onto the rollator handle tubing. Feed the cable through the leg bracket, then feed the spring onto the cable as shown. Next support the cable and slide the brass cable stop onto the end of the cable while trying to avoid bending any cable strands.

Once the brass cable stop is flush with the spring, tighten the socket head cap screw.

5. Support the brake shoe, lift the spring, and compress it enough so the shoe will fit between the spring and the brass cable stop. The cable stop must go completely below the brake shoe.

Notice, at this point the brake is not tensioned properly and the brake shoe is far from the tire.



6. Make sure the brake lever is NOT in the locked position. While supporting the cable with the pliers, loosen the locking screw on the brass cable stop.

Next lift the cable stop up while pulling down on the cable with pliers. There should be about 2 millimeters of space between the brake shoe and the tire. (Thickness of two pennies stacked). Confirm the brake is functional and adjust if needed.

7. You must use adhesive to reinstall the grip. Spray is recommended but epoxy, gorilla glue and rubber cement are suitable adhesives.

Apply a moderate amount of spray adhesive (or a few drops of a recommended alternative glue) to handle tubing.

8. Position the handles at the same height and secure with knob and bolt. Slide the grip over the adhesive completely to the brake assembly. Before allowing the adhesive to set, grasp the handles as if you are operating the rollator and twist to ensure they are positioned correctly and feel the same on both sides.

9. After following the recommended instructions to apply the adhesive, allow adhesive to set before using the rollator.

\* It is recommended that you seek the assistance of a local provider in replacing your rollator handbrakes if and when possible.

