



The Zenith II™ Full-Electric Bed

Carefully designed to meet important clinical requirements of Long Term Care. From a low of 7.9 inches to a high position of 30 inches, the Zenith II bed facilitates resident safety and ensures caregiver ease of operation. Driven by quiet, powerful DC motors, its innovative patent-pending design is engineered to keep the bed completely level at all times. With its improved strength, faster hi-lo travel, and superior features, it's the only bed you'll ever need, and will fit reasonably in your budget.

Ask about our **FAST** Hi-Lo Travel Options!

- Conforms to UL 60601-1 and IEC601-2-38; these are the highest Safety Standards for Medical Beds in Long Term Care, including:
 - No pinch points within 200mm (approx. 8") of the bed perimeter.
 - Rolling clearance of 20mm (approx. 0.8") over thresholds.
 - Spash-proof electrical connections.
 - Load capacity and tip-stability testing.
- Strong and swift vertical travel with 500 lb safe working load capacity.
- Locks or rolls with 360° mobility at any height, even the lowest.
- Rigid tubular cross-members beneath grid or pan mattress deck for maximum support.
- Assist Devices remove w/o tools, and meet FDA Entrapment Guidelines.

76" LENGTH	80" LENGTH	Mattress Deck Surface	Bed Type	Units
ZZ6074	ZZ6174	Grid	Full-Electric w/Z-lok ²	1 ea
ZZ6014	ZZ6114	Pan	Full-Electric w/Z-lok ²	1 ea
ZG6074	ZG6174	Grid	Full-Elect. w/Standard Caster	1 ea
ZG6014	ZG6114	Pan	Full-Elect. w/Standard Caster	1 ea
ZG6074S	ZG6174S	Grid	Full-Elect. w/Steel Caster	1 ea
ZG6014S	ZG6114S	Pan	Full-Elect. w/Steel Caster	1 ea

Bed, head and foot boards, and rails sold separately.





Zenith™ II Full-Electric Bed shown with Z-Lok Pedals, Grid Mattress Deck, Stock Windsor Mahogany Center Crown Head Board, Stock Windsor Mahogany Center Crown T-Mold Foot Board with Embedded Staff Control, and Head Half Rails.



Optional Steel Caster



Pendant Holster

Options and Accessories

Head and Foot Boards*	
Z53500	Rectangular, Solar Oak Finish
Z48500	Rectangular, Medium Oak Finish
Z80500	Rectangular, Montana Walnut Finish
ZL788053	Rectangular, Solar Oak Finish, Embedded Staff Control
ZL788048	Rectangular, Medium Oak Finish, Embedded Staff Control
ZL788080	Rectangular, Montana Walnut Finish, Embedded Staff Control
Z29500CCR	Center Crown, Windsor Mahogany Finish
Z02500CCR	Center Crown, Figured Mahogany Finish
ZL788029CCR	Center Crown, Windsor Mahogany Finish, Embedded Staff Control
ZL788002CCR	Center Crown, Figured Mahogany Finish, Embedded Staff Control

Staff Control Electronics Pkg w/o Boards	
ZL788000	For use w/custom boards

* Other finishes available upon request.
All boards above have T-mold edges.

Assist Devices and Assist Bars		
76" LENGTH	80" LENGTH	Descriptions
ZA78400	ZA78400	Half Length Assist Head End (1Pr)
ZA79400	ZA78500	Half Length Assist Foot End (1Pr)
Head End	Foot End	Assist Bars
ZA78600	ZA78650	Assist Bar, 76" & 80" (ea)

Wall Saver Bumpers	
ZA78700	Widesaver Tubular Bumper-Head End
ZA78300	Side Wall Saver Bumper (Set of 2)

Accessories	
ZA78900	Economy Lockout Box
ZA71900	4" Mattress Deck Extension
ZA79000	Trapeze Adapter (For use with ZA78100)
ZA78100	Trapeze Bar-Offset (must use with ZA79000)
ZA76800	Traction Adapter Kit
999-0711-003	Back-Up Battery Pack
999-0711-002	Battery Charger
554-0769-000	Corner Mattress Retainer Kit
554-0190-000	Board Mounted I.V. Sockets

