

Bunn BA 400A Silent Air Compressor

- Provides air source for majority of aerosol generators and croup tents.
- Designed for efficient, continuous, high-pressure performance.
- Unique design of WOB-L® piston compressor provides reliable and quiet 48 dBA operation.
- Vinyl finish cabinet for quick clean-up.
- Easy to read gauges.
- Locking regulator knob for consistent delivery.
- Unique internal mounting bracket for easy set-up and transport.
- Compact size permits storage on stands designed for aerosol tent applications.
- Designed to meet the requirements of HCPCS code: E0565.
- One year limited warranty.

JB0120-010 1 ea



Technical Specifications

Width	14"
Depth	12"
Height	12"
Weight	20 lb
Power Requirements	3.5 A @ 115 VAC, 60 Hz
Compressor Pressure Range	0 – 100 psi
Compressor Flow Range	0 – 35 lpm



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com