



Matrix Plus™

The Matrix Plus comes with many standard features, including sleep deck height adjustment from 8.75" to 25.25" and a single motor for Hi-Lo function that ensures the bed is always level. It has a dynamic operating capacity of 450 lb and meets the highest UL safety standards for medical beds, UL60601 and EN60601. CAN/CSA C22.2 No. 601-1-M90 approved. It has strong tubular panel mounting and built-in stabilizer feet that ensure the bed will automatically lock in the lowest position. The bed features individual locking casters. The bed system is patent pending.

Specifications:

- Sleep deck dimensions: 35" width, and 76" or 80" sleep deck lengths.
- Maximum head elevation of 70 degrees.
- Sleep surface constructed of strong steel wire grid; no plastic or springs are used.
- Rigid construction utilizing heavy tubular steel to maximize strength and stability of the sleep surface in all positions.
- All weldments and components are electrostatically powder-coated or plated.
- 115-Volt standard. 230 Volt available
- Full-Electric operation with three ultra whisper quiet, rapid-moving DC motors.
- Standard with handheld 6-function control.

Warranty:

- Electronics: 3 Years
- All other components: 15 Years

76" Length	80" Length	Sleep Surface	Bed Type	Units
MP5074	MP5174	Grid Deck	Full-Electric	1 ea
MP5014	MP5114	Pan Deck	Full-Electric	1 ea
MM5074	MM5174	Grid Deck	Full-Elect. w/ M-Lok	1 ea
MM5014	MM5114	Pan Deck	Full-Elect. w/ M-Lok	1 ea

* Bed, head and foot boards, and siderails sold separately.





Matrix Plus Bed shown with Grid Deck and Custom Windsor Mahogany Center Crown Wood Cap Headboard, Stock Windsor Mahogany Center Crown T-Mold Footboard with Embedded Staff Control, and Stabilizer Legs.

Options and Accessories

Head and Foot Boards*	
Z46500	Rectangular, Light Oak Finish
Z48500	Rectangular, Medium Oak Finish
Z80500	Rectangular, Dark Walnut Finish
Z46550	Rectangular, Light Oak Finish, Embedded Staff Control**
Z48550	Rectangular, Medium Oak Finish, Embedded Staff Control**
Z80550	Rectangular, Dark Walnut Finish, Embedded Staff Control**
Z29500CCR	Center Crown, Windsor Mahogany Finish
Z29550CCR	Center Crown, Windsor Mahogany Finish, Embedded Staff Control**
Z02500CCR	Center Crown, Figured Mahogany Finish
Z02550CCR	Center Crown, Figured Mahogany Finish, Embedded Staff Control**

* Other finishes available upon request.

** Order in conjunction with ZL74700 Embedded Staff Control.

*** All boards above have T-mold edges.

Side Rails and Assist Bars	
ZA7020	Half Head Rails (Set) or Foot Siderails
AZ7350	Assist Bar (1 Each)

* Foot rails for 80" bed only.

Factory-Installed Options	
ZL71510	Single-Function Lock-Out
ZL71530	Three-Function Lock-Out
ZL74700	Embedded Staff Control
999-0711-005	Grounded Hospital-Grade Plug

Accessories	
999-0711-003	Back-Up Battery Pack
999-0711-002	Back-Up Battery Charger
554-0143-000	Trapeze Bar (Patient Helper Bar)
554-0190-000	Panel-Mounted I.V. Holder
554-0769-000	Corner Mattress Retainer Kit
AM75500	Tubular Wall Saver
554-0720-000	Trapeze Adaptor
554-0743-000	Pendant Holster
554-0721-000	Quick Release Head or Foot Board Mounting Kit
554-0768-000	Traction Adapter Kit
554-0767-000	Edema Foot Ratchet
554-0719-000	4" Sleep Deck Extension

