

## Lightweight Rigid Walker

- Non-folding design provides extra rigidity for use where space is not an issue.
- Angled legs provide maximum patient support in all directions and help resist tipping.
- Angled side brace helps support patients as they stand up and sit down.
- Fits easily around standard toilet fixtures.
- Durable 1" aluminum tubing provides maximum strength while remaining lightweight.
- Adjustable to fit a broad range of user heights.
- Accommodates all leg attachments.
- Designed to meet the requirements of HCPCS code: E0130.
- Limited lifetime warranty.

**1925A** 4/cs



]	Item	Width Inside Hand Grips	Width Inside Rear Legs	Overall Width	Depth Front to Rear*	Minimum Height	Maximum Height	Approx. User Height	Weight Capacity	Units Per Case
1	925A	17"	19 ¾"	22 ½"	20 ½"	31"	38"	5' 3" - 6' 2"	400 lb	4

