

Halogen Exam Lamp

- Flared, double-wall, dual cool ventilated reflector and shade are of heavy gauge anodized aluminum.
- 50 watt bright, focused, white halogen light
- Designed with swivel assembly for complete adjustability.
- One-hand, slip proof lock. Weighted, chrome-plated steel base.
- 115/120V AC transformer 12V 50W.
- Height adjustable from 34 ½" – 53 ½".

1696	Grafcó Halogen Exam Lamp	1 ea
1696M	Grafcó Halogen Exam Lamp, Mobile Base	1 ea
1696-50HAL	Replacement Halogen Bulb For 1696, 1696M	1 ea



1696M



1696



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0700147RevA07