

HybridLX Rollator Transport Chair

The HybridLX is an innovative rollator that combines the function of a wheeled walker and a transport chair into one lightweight design. Keeping the patient's desire for independence and comfort in mind, the HybridLX is designed with secure handgrips, a comfortable and wide seat and backrest, and the effortless design to transition from a rollator to a transport chair.

Other Features Include:

- Lightweight frame – durable, aluminum frame weighs only 22 lbs (this includes footrests) making it easier for users to move about freely
- Hybrid concept – designed to be used primarily as a rollator while easily transforming into a transport chair with removable footrests (adjustable in length) when user needs a short rest
- Ultra-Compact design – easy to assemble without the need for tools and folds effortlessly for convenient transport and storage
- Easy-to-use handbrakes – allow for short and longer period braking when used as rollator and transport chair (user must be accompanied when using as a transport chair)
- Flexible, soft polyester backrest and 18" wide mesh seat with pull up strap for easy folding provides comfort for user
- Large 8" wheels for easy rolling indoor and outdoor on rough surfaces
- Complete with storage pouch for personal items
- Available in attractive, contemporary colors: Titanium and Majestic Blue
- Designed to meet the requirements of HCPCS Codes: E0143 + E0156.
- Warranty: 3 years on Frame, 1 year on other components

LX1000T	Titanium	1 ea
LX1000B	Majestic Blue	1 ea
LX1000-CB	Carry Bag	1 ea



Specifications

Width Inside Handgrips	18.5" (470mm)
Overall Width	26" (660mm)
Overall Depth (with footrests)	40" (1025mm)
Depth Front to Rear	28.5" (730mm)
Seat Width	18" (460mm)
Seat Depth	11" (280mm)
Seat to Floor	21.25" (540mm)
Wheels	8" Black
Minimum Handle Height	30" (770mm)
Maximum Handle Height	37" (945mm)
Approx. User Height	5'4" - 6'2"
Weight Capacity	300 lbs
Product Weight	22 lbs (9.88kg) with footrest



Graham-Field
2935 Northeast Parkway
Atlanta, Georgia 30360
© 2008 GF Health Products, Inc.

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0800142RevA08