

Multi-Cuff-5 Blood Pressure Kit

- Compact orange multiple blood pressure cuff kit including nylon infant, child, adult, large adult, and thigh blood pressure cuffs
- Durable, aluminum manometer with large, easy-to-read 65mm dial, 1-tube inflation bladder, integral latex-free, ribbed inflation bulb with dust filter, and deluxe push button air release valve
- The rotating bulb and valve creates the ability for ambidextrous use
- Includes extra storage compartment with hook and loop type closure

206 Orange 1 ea

