

Traveler® SE

Everest & Jennings introduces the latest in a long tradition of quality wheelchairs. The Traveler SE is the newest addition to the Traveler line; a brand that's well known and respected for performance and durability.

The Traveler SE is the standard wheelchair of choice. Designed with dual rear axles and adjustable casters, this wheelchair can convert from a K0001 Medicare code to a K0002 Medicare code. The Traveler SE has industry compatible front rigging and has vinyl upholstery constructed with a strong inner lining for added support.

Standard Features

- Black vinyl upholstery
- Chart pocket on back
- "Tool Free" adjustable front rigging
- Hammer tone paint provides a long lasting and durable finish
- Industry compatible front rigging
- Two position axle and caster fork provides 2" seat height adjustment
- 250 lb. weight capacity
- Meets Medicare code K0001 and K0002

Warranty

- Limited 1 year on side frames and crossbars



Traveler® SE		Dimensions and Weights								
		Seat			Back	Overall				Wt. *
3E000000	Fixed Full Arm, Fixed Footrest	18"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	26 1/2"	12"	40 lb
3E010100	Fixed Full Arm, Swingaway Footrest	18"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	26 1/2"	12"	40 lb
3E010110	Fixed Full Arm, Elevating Legrest	18"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	26 1/2"	12"	40 lb
3E010120	Detachable Desk Arm, Swingaway Footrest	18"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	26 1/2"	12"	40 lb
3E010130	Detachable Desk Arm, Elevating Legrest	18"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	26 1/2"	12"	40 lb
3E010140	Detachable Full Arm, Swingaway Footrest	18"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	26 1/2"	12"	40 lb
3E010150	Detachable Full Arm, Elevating Legrest	18"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	26 1/2"	12"	40 lb
3E010200	Fixed Full Arm, Swingaway Footrest	16"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	24 1/2"	12"	40 lb
3E010210	Fixed Full Arm, Elevating Legrest	16"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	24 1/2"	12"	40 lb
3E010220	Detachable Desk Arm, Swingaway Footrest	16"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	24 1/2"	12"	40 lb
3E010230	Detachable Desk Arm, Elevating Legrest	16"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	24 1/2"	12"	40 lb
3E010240	Detachable Full Arm, Swingaway Footrest	16"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	24 1/2"	12"	40 lb
3E010250	Detachable Full Arm, Elevating Legrest	16"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	24 1/2"	12"	40 lb
3E010320	Detachable Desk Arm, Swingaway Footrest	20"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	28 1/2"	12"	40 lb
3E010330	Detachable Desk Arm, Elevating Legrest	20"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	28 1/2"	12"	40 lb
3E010340	Detachable Full Arm, Swingaway Footrest	20"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	28 1/2"	12"	40 lb
3E010350	Detachable Full Arm, Elevating Legrest	20"	16"	18" or 20"	15 1/2"	31 1/2"	35 1/2"	28 1/2"	12"	40 lb

NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4". *Without Front Rigging.

