

## Bariatric Recliner

The Lumex Electric Bariatric Recliner is the perfect solution for larger patients who require a wider chair. Composed of a heavy-gauge steel frame, this chair has a maximum weight capacity of 700 lb! A four-function hand pendant activates the motor allowing caregivers to safely recline and raise patients without strain.

- Quick-release trigger mechanism on push handle allows chair to swiftly move into Trendelenburg position.
- Fully equipped with rechargeable battery pack for full operation away from a wall outlet.
- Four swivel casters facilitate transport and maneuverability in confined areas.
- Durable, high impact side tables sit flush with chair arm to conveniently hold prep materials. Easy to activate by caregiver.
- Legrest and backrest raise and lower independently, giving the operator complete control over positioning.
- Comes complete with attaching head pillow for extra support and added comfort.
- Limited lifetime warranty on frame.
- One year warranty on electronics.
- Three years warranty on all other components.


### Specs

Overall Height.....	55"
Overall Width.....	38"
Overall Depth.....	59"
Width Between Arms.....	25.25"
Arm Height From Floor.....	32"
Arm Height From Seat.....	8"
Seat Height.....	25.5"
Seat Width.....	25.75"
Seat Depth.....	24"
Depth Fully Reclined In Trendelenburg .....	71"
Width With Both Side Tables in Upright Position .....	58"
Shipping Weight .....	308 lb
Max. Weight Capacity.....	700 lb

### Available Colors:

<b>FR588W427</b>	Blue Ridge	1 ea
<b>FR588W857</b>	Jade	1 ea
<b>FR588W863</b>	Rosewood	1 ea

(Meets California Technical Bulletin 133  
Flammability Standards)



Blue Ridge (427)



Jade (857)



Rosewood (863)



Lumex Clinical Care Recliners are equipped with LiquiCell and AEGIS Microbe Shield® technologies.

LiquiCell is an ultra-thin, liquid-filled interface specifically designed and proven to significantly reduce shear stress, aiding in the prevention of pressure sores and general discomfort.

AEGIS Microbe Shield works by chemically bonding itself to the treated surface and rupturing the membranes of any bacteria that it comes in contact with, thereby inhibiting microbial growth.

- Upholstery color may vary slightly from brochure.
- Lumex Ortho-Biotic Seating with an FR prefix designates that the seating meets CA133 Fire Code specifications.
- Please contact Graham-Field Customer Service at 800-347-5678 for details on specialty colors and COM orders.



GRAHAM-FIELD

GF0900052RevA09

**www.grahamfield.com**

© 2009 **GF Health Products, Inc.**  
2935 Northeast Parkway  
Atlanta, Georgia 30360