

Drop-Arm Recliners (Model 3302)

- Pillow back, three-positions, fire-rated upholstery.
- Combines the therapeutic benefits of position change with a more residential appearance. Designed to aid in circulation and to help ease discomfort from respiratory, heart and other ailments.
- Unique drop-arm mechanism facilitates lateral transfers. Easy to activate by caregiver.
- Ideal for knee and hip surgery patients.
- Three inch ball bearing casters ensure smooth transport; foot actuated swivel locking rear casters for a more secure lock.
- Unitized welded steel frame for durability. Flame and stain resistant vinyl upholstery available in a variety of colors.
- Weight Bearing Capacity: 250 Pounds.

• Lumex Wrap-Around Warranty

Welded Frame: Limited Lifetime

Upholstery: Limited One Year

All Other Components: Limited Three Year

Available Colors:

FR3302427	Blue Ridge	1 ea
FR3302851	Doe Skin	1 ea
FR3302857	Jade	1 ea
FR3302863	Rosewood	1 ea

*(Meets California Technical Bulletin 133
Flammability Standards)*



- Upholstery color may vary slightly from brochure.
- Lumex Ortho-Biotic Seating with an FR prefix designates that the seating meets CA133 Fire Code specifications.
- Please contact Graham-Field Customer Service at 800-347-5678 for details on specialty colors.



GRAHAM-FIELD

GF0900053RevA09

www.grahamfield.com

© 2009 **GF Health Products, Inc.**
2935 Northeast Parkway
Atlanta, Georgia 30360