

## Economy Overbed Table, Non-Tilt

- Easy to assemble
- Height infinitely adjusts 14"
- Swivel casters provide smooth-rolling action. Two locking and two non-locking casters
- Walnut, wood-grain laminated 31" top on ¾" particle-board core with plastic edge molding
- Chrome plated "H" base and column for aesthetics and durability
- Designed to meet the requirements of HCPCS code: E0274
- Supports up to 30 lbs. of dispersed weight.
- Warranty : One Year

**GF8900-1** 1 ea



### Specifications

Table Height .....	31.25" - 42.25"
Table Width .....	31.5"
Table Depth .....	15"
Base Width .....	26.2"
Base Depth .....	16.2"
Weight .....	22 lbs
Dispersed Weight Capacity .....	30 lbs



**GRAHAM-FIELD**

GF0700105RevB09

**[www.grahamfield.com](http://www.grahamfield.com)**

© 2009 **GF Health Products, Inc.**  
2935 Northeast Parkway  
Atlanta, Georgia 30360