

Preferred Care® Recliner Series

The Lumex 565 Patented recline mechanism allows for infinite positioning of recline. Reduced armrest and seat height enables the recliner to fit under dining/activity tables and allows resident to "scoot" and more easily reposition chair by themselves. Two front swivel casters and two rear locking casters enable the Preferred Care Recliner to be easily maneuvered in small spaces and moved down hallways.

Standard Features

• Infinite recliner back patented mechanism provides infinite positioning of the back, ranging from 95 degrees to 145 degrees for increased resident comfort and pressure relief. The reclining back can be locked in any position by the caregiver.

• A unique suspension system self-adjusts to each person, cradling the resident from head to heel and virtually eliminating the risk of "bottoming out".

• Extra-wide padded armrests and side panels maximize comfort.

 Durable epoxy coated welded steel frame is easy to clean and maintain.

Quick-release upholstery allows for easy cleaning, removal and replacements. Vinyl upholstered back is sulfide stain, mildew and oil resistant.

• A variety of accessories are available to customize the chair for specific individual needs.

• Lumex Wrap-Around Warranty Welded Frame: Limited Lifetime **Upholstery:** Limited One Year

All Other Components: Limited Three Year

Available Colors:

565G401	Vintage Gold	1 ea
565G409	Warm Taupe	1 ea
565G427	Blue Ridge	1 ea
565G432	Imperial Blue	1 ea
565G851	Doe Skin	1 ea
565G857	Jade	1 ea
565G863	Rosewood	1 ea

(Meets California Technical Bulletin 133

Flammability Standards)

FR565G427	Blue Ridge	1 ea
FR565G857	Jade	1 ea
FR565G863	Rosewood	1 ea

Doe Skin (851)



Tray Table and bolsters sold separately and are not included with recliner.



Jade (857)

Rosewood (863)

- Upholstery color may vary slightly from brochure.
- Lumex Ortho-Biotic Seating with an FR prefix designates that the seating meets CA133 Fire Code specifications.
- Please contact Graham-Field Customer Service at 800-347-5678 for details on specialty colors and COM orders.

