



# INTRO BED ORDER FORM

Form# 999-0816-185  
REV A 02/10

## Model 705

To order **NEW** beds, please contact a GF Health Products Sales Representative at 1-800-347-5678.

### INTRO 705 BEDS AVAILABLE

Order #	Description	Qty Ordered	Order Date
<b>NT Series with Hi/Lo Motor, 3" Twin Locking Casters, and Wheels</b>			
NT7074	76" Grid Deck — Full electric with Head, Foot, and Hi/Lo Motor		
NT7174	80" Grid Deck — Full electric with with Head, Foot, and Hi/Lo Motor		
NT7014	76" Pan Deck — Full electric with with Head, Foot, and Hi/Lo Motor		
NT7114	80" Pan Deck — Full electric with with Head, Foot, and Hi/Lo Motor		
<b>NT Series with Fixed Height and 3" Twin Locking Casters</b>			
NT7077	76" Grid Deck — Full electric with Head and Foot Motor		
NT7177	80" Grid Deck — Full electric with Head and Foot Motor		
NT7017	76" Pan Deck — Full electric with Head and Foot Motor		
NT7017	80" Pan Deck — Full electric with Head and Foot Motor		

**NOTE: Footboard with embedded Staff/Nurse Control is optional and is not part of the standard Intro Bed.**

### TABLE OF CONTENTS:

For new or replacement bed service parts please refer to the following pages.

<b>Main Bed Section</b> (Decking, u-levers, radial arm tube, caster carriage, casters, etc.)	Page 2
<b>Boards With Staff Control Assemblies</b> (Hd/Ft Boards, Staff Control, Cables, Shroud, etc.)	Page 2
<b>Bed End Components</b> (Standard boards, embedded footboards, mounting kits, etc).	Page 3
<b>Electronic Components</b> (Cables, motors, control boxes, hand controllers/pendants, etc)	Page 3
<b>Wallsavers and Mattress Retainers</b>	Page 4
<b>Assist Devices</b> (Half rails, foot mounting bracket, assist bars, trapeze adaptors, trapeze supports)	Page 4
<b>Accessories</b> (4" extension kits, wide kits, I.V. mount, inclinometer)	Page 4
<b>Economy Lock Out &amp; Battery Backup Options</b>	Page 5
<b>Bed Hardware</b>	Page 5
<b>Control Box and Motor Mounting Hardware</b>	Page 6
<b>Miscellaneous Hardware and Literature</b>	Page 6
<b>Standard Intro 705 Bed Diagram</b>	Page 7



# INTRO BED SERVICE PART ORDER FORM

**Model 705**

Form# 999-0816-185  
REV A 02/10

The following list of service parts for the Intro 705 Bed have been provided for your convenience. Please keep this list in a secure location and use as a reference when ordering new or replacement parts. A general bed diagram has been provided on the last page.

**To order new or replacement SERVICE PARTS, please contact a Basic American Medical Products Customer Service Representative at 1-800-365-2338 to confirm you are ordering the right part.**

MAIN BED SECTIONS		Please see back diagram (page 7) for general part identification.			
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
DECKING					
Same as Part#	999-0816-900G	Main Frame Weldment (Reference-NOT A SERVICE PART)	1		
Same as Part#	999-0779-920G	Standard Head/Back Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0779-930G	Standard Foot Grid Deck Weldment (80")	1		
Same as Part#	999-0775-936G	Optional Foot Grid Deck Weldment (76")	(1)		
Same as Part#	999-0779-940G	Standard Knee Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0796-982G	Optional Head/Back Pan Deck Weldment (76" & 80")	(1)		
Same as Part#	999-0796-983G	Optional Foot Pan Deck Weldment (80")	(1)		
Same as Part#	999-0796-986G	Optional Foot Pan Deck Weldment (76")	(1)		
Same as Part#	999-0796-984G	Optional Knee Pan Deck Weldment (76" & 80")	(1)		
Same as Part#	999-0777-022G	Long Main Frame Mounting Bracket	4		
Same as Part#	999-0807-036G	Leveler Strap Bracket	2		
HEAD & FOOT U-LEVER/CASTER BASES & DRAG LINK					
Same as Part#	999-0778-903G	Head U-Lever/Caster Base Weldment	1		
Same as Part#	999-0778-904G	Foot U-Lever/Caster Base Weldment	1		
Same as Part#	999-0775-906G	Drag Link Weldment	1		
CASTERS & WHEELS (if ordered separately)					
Same as Part#	100-6330-061	Single 3" Non-Locking Casters for Standard Carriage	2		
Same as Part#	100-6330-062	Single 3" Locking Casters for Standard Carriage	2		
Same as Part#	100-6325-003	Single 2.50" Dia. Rubber Wheel	4		
CAPPING (if ordered separately)					
Same as Part#	999-0775-001	Half Moon Caps (Head & Foot Caster Bases/U-Levers)	12		
Same as Part#	100-4700-015	1.00" Square Tube End Cap (Main Frame)	4		
Same as Part#	100-4700-017	1.00" x 2.00" Rect. Tube End Cap (2 Main Frame & 4 Carriage)	2		

HD/FT BOARDS WITH EMBEDDED STAFF CONTROL ASSEMBLIES					
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
ZL788053	-	Solar Oak Hd/Ft. Boards with Staff Control	1/1/1		
ZL788048	-	Med. Oak Hd/Ft. Boards with Staff Control	1/1/1		
ZL788080	-	Montana Walnut Hd/Ft. Boards with Staff Control	1/1/1		
ZL788029CCR	-	Windsor Mahogany Hd/Ft. Boards with Staff Control	1/1/1		
ZL788002CCR	-	Figured Mahogany Hd/Ft. Boards with Staff Control	1/1/1		
ZL788004CCR	-	Wild Cherry Hd/Ft. Boards with Staff Control	1/1/1		

**BED END COMPONENTS**

Please see back diagram (page 7) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
<b>STANDARD HEAD- &amp; FOOTBOARDS WITHOUT CUTOUT:</b>					
Z53500	-	Solar Oak Hd/Ft Boards without staff control cutout	1/1		
Z48500	-	Med Oak Hd/Ft Boards without staff control cutout	1/1		
Z80500	-	Walnut Oak Hd/Ft Boards without staff control cutout	1/1		
Z29500CCR	-	Wind Mah CC Hd/Ft Boards without staff control cutout	1/1		
Z04500CCR	-	Wild Cher CC Hd/Ft Boards without staff control cutout	1/1		
Same as Part#	100-7025-004	1/4-20 UNC x .60 Long Threaded Inserts (pre-inserted into footboard and headboard)	8 (4 per bd)		
<b>STANDARD HEAD- &amp; FOOTBOARDS WITH STAFF CONTROL CUTOUT:</b>					
Z53650	-	Solar Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z48650	-	Med Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z80650	-	Walnut Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z29650CCR	-	Wind Mah CC Hd/Ft Boards cutout for emb. staff control	1/1		
Z04650CCR	-	Wild Cher CC Hd/Ft Boards cutout for emb. staff control	1/1		
Same as Part#	100-7025-004	1/4-20 UNC x .60 Long Threaded Inserts (pre-inserted into footboard and headboard)	8 (4 per bd)		
<b>STANDARD MOUNTING KIT 999-0775-000 INCLUDES:</b>					
Same as Part#	999-0746-035G	Angled Mounting Tubes (for mounting boards to frame)	4		
Same as Part#	100-4700-015	1" Square Tube End Caps (pre-inserted in mounting tubes)	4		
Same as Part#	100-5425-081	1/4-20 UNC x 1.50" Phillips Pan Head Machine Screws	8		
Same as Part#	100-7931-007	5/16 x 1.50" Clevis Pin with hole (mounting tubes to frame)	4		
Same as Part#	100-2001-006	1-9/16 Hair Pin Clip (mounting tubes to frame)	4		
Same as Part#	999-0775-195	Assembly Instruction Sheet	1		

**ELECTRICAL COMPONENTS**

Please see back diagram (page 7) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	999-0775-052	Head Actuator/Motor	1		
Same as Part#	999-0778-204	Head Actuator/Motor Cable (8")	1		
Same as Part#	999-0775-053	Foot/Knee Actuator/Motor	1		
Same as Part#	999-0775-202	Foot/Knee Actuator/Motor Cable (60")	1		
Same as Part#	999-0778-051	Standard Hi/Lo Actuator/Motors	1		
Same as Part#	999-0775-201	Hi/Lo Actuator/Motor Cable (28")	1		
Same as Part#	999-0779-300	J-Box/Control Box	1		
Same as Part#	999-0775-208	J-Box/Control Box 3-Prong Power Cable with Ground	1		
Same as Part#	999-0775-305	Standard Pendant/Hand Controller	1		
Same as Part#	999-0791-000	Pendant Holster	1		
Same as Part#	999-0779-200	Optional Control Box/Staff Control Extension Cable (24")	1		
-	999-0775-901	Optional Staff Control Assembly (with 2 cables, mounting bezel, and control panel)	1 Assm.		
Same as Part#	999-0775-006G	Optional Wide Shroud (covers staff control cables on footboard)	1		
Same as Part#	100-1100-024	4" Cable Tie (staff control cables to footboard mounting tube)	1		
Same as Part#	100-5506-002	#6 x .50" Phillips Truss Head Screws (shroud & control)	6		
ZL788000	-	Electronics Kit for embedded staff control - includes Control Panel, Bezel, Attached Cables, Shroud, Instruction Sheet, attaching screws, and 4" cable tie	1 kit		
Same as Part#	100-1100-018	Cable Ties for all cables to frame	22		

OTHER COMPONENTS			Please see back diagram (page 7) for general part identification.		
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
ZA79800	-	Optional Intro Wallsaver Kit (Includes standard size wallsaver, 999-0798-001, mounting hardware & assembly instructions)	1 kit		
ZA79820	-	Trapeze Intro Wallsaver Kit (Includes trapeze size wallsaver, 999-0798-200, mounting hardware & assembly instructions)	1 kit		
999-0814-000	-	Corner Wire Mattress Retainer Kit (set of two retainers)	1 kit		
Same as Part#	999-0814-001	Single Wire Mattress Retainer	2		

ASSIST DEVICES		NOTE: Rail "Kits" include both Quick Release Lanyard and Locking Bolt/Nut mounting hardware			
Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA78400 (both rails & hardware)	-	Left - Head Half Side Rails (76" & 80") Right - Head Half Side Rails (76" & 80")	1/1/4 for head end		
ZA78500 (includes both rails & two foot bracket kits)	-	Left - Foot Half Side Rails (80") with hardware Right - Foot Half Side Rails (80") with hardware Foot Bracket Kits with hardware	1/1/4 for 80" foot end		
ZA79400 (includes both rails & two foot bracket kits)	-	Left - Foot Half Side Rails (76") with hardware Right - Foot Half Side Rails (76") with hardware Foot Bracket Kits with hardware	1/1/4 for 76" foot end		
	999-0784-907G	Rt. Head (76"/80") or Lft. Foot (80") Half Rail only	1		
	999-0784-908G	Lft. Head (76"/80") or Rt. Foot (80") Half Rail only	1		
	999-0794-907G	Left Foot Half Rail for 76" only	1		
	999-0794-908G	Right Foot Half Rail for 76" only	1		
-	999-0786-995	Foot Bracket Hardware Bag with Bracket	1		
-	999-0775-025G	Foot Bracket only	1		
ZA80300 (includes hardware)		One Assist Bar with Mounting Hardware (Includes Quick Release Lanyard and Bolt/Nut Hardware)	1/2 for head end		
ZA80350 (includes 1 foot bracket)	-	One Assist Bar with mounting hardware, and Foot Bracket Kit with mounting hardware	1/2 for foot end		
-	999-0803-902G	One Assist Bar only (without hardware)	1		
-	999-0784-100	Quick Release Lanyard only	1		
ZA79000	-	Trapeze Support Adaptor Kit (Includes Trapeze Support and attaching hardware)	1 for head end		
ZA78100	-	Trapeze Bar Offset (must use with ZA79000)	1		

ACCESSORIES					
Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA71900	-	4" Standard Extension Foot Pan Kit	1 kit		
ZA82100		4" Extension Foot Pan Kit for Wide Deck	1 kit		
ZA76800	-	Traction Adaptor Kit	1		
ZA19000	554-0190-000	Brd. Mounted I.V. Sockets (includes bracket, cap & hdwr)	1		
WIDE KITS FOR GEN 7 INTRO BEDS					
ZA811000	-	42" Wide Kit for Rail Mount - No Boards			
WIDE KITS FOR GEN 7 INTRO BEDS WITH BOARDS & STAFF CONTROL ASSEMBLIES					
ZA811053	-	42" Wide Kit for Rail Mount - Solar Oak Bds w/Staff Control			
ZA811048	-	42" Wide Kit for Rail Mount - Med Oak Bds w/Staff Control			
ZA811080	-	42" Wide Kit for Rail Mount - Walnut Oak Bds w/Staff Control			
ZA811004CCR	-	42" Wide Kit for Rail Mount - CherCC Oak Bds w/Staff Control			
ZA811029CCR	-	42" Wide Kit for Rail Mount - WinCC Oak Bds w/Staff Control			

<b>ECONOMY LOCK OUT AND BACKUP BATTERY CHARGER</b>					
<b>Order #</b>	<b>Part #</b>	<b>Description</b>	<b>Qty On Bed</b>	<b>Qty Ordered</b>	<b>Order Date</b>
ZA78900	-	Economy Lockout Box	1		
Same as Part#	999-0711-002	Battery Charger/Transformer (order with battery pack)	1		
Same as Part#	999-0711-003	Backup Battery Pack (order with charger/transformer)	1		

<b>BED HARDWARE</b>					
<b>Purchase #</b>	<b>Part #</b>	<b>Description</b>	<b>Qty On Bed</b>	<b>Qty Ordered</b>	<b>Order Date</b>
<b>DECK AND LEVELER STRAPS TO FRAME MOUNTING HARDWARE</b>					
Same as Part#	100-5431-002	5/16-18 x .375" Hex Washer Head Machine Screw Secures Head & Knee Decks to Pan Seat and Knee to Foot Deck, and Leveler Straps to Foot Deck Brackets	8		
Same as Part#	999-0032-010	Short Self-Wrenching Nut Secures Head & Knee Decks to Pan Seat and Knee to Foot Deck	6		
Same as Part#	100-7950-004	.50" x 2.00" Clevis Pin with hole Secures Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-7518-001	7/16-.50" Circle Cotter Secures Clevis Pin for Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-5231-005	5/16-18 x 1.75" Carriage Bolt Secures Leveler Straps to Main Frame Weldment Rails	2		
Same as Part#	100-6731-006	5/16-18 Thin Hex Head Nylon Lock Nut Secures Bolts for Leveler Straps to Main Frame Rails	2		
Same as Part#	806-0277-008	Long Self-Wrenching Nut Secures Leveler Straps to Foot Deck Brackets	2		
<b>LONG MOUNTING BRACKETS &amp; HEAD MOTOR BRACKET TO MAIN FRAME</b>					
Same as Part#	100-5231-005	5/16-18 x 1.75" Carriage Bolt Secures Mounting Brackets & Motor Bracket to Main Frame Rails	9		
Same as Part#	100-6731-006	5/16-18 Thin Hex Head Nylon Lock Nut Secures Mounting Brackets & Motor Bracket to Main Frame Rails	9		
<b>DRAG LINK TO MAIN FRAME AND RADIAL ARM MOUNTING HARDWARE</b>					
Same as Part#	999-0775-450	12mm x 50mm Clevis Pin w/hole Secures Drag Link Brackets to Head- & Foot Caster Base Brackets	2		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Drag Link Clevis Pins to Head- & Foot Caster Base Brackets	2		
<b>HEAD &amp; FOOT CASTER BASE/U-LEVER WELDMENTS TO MAIN FRAME MOUNTING BRACKETS</b>					
Same as Part#	999-0775-450	12mm x 40mm Clevis Pin w/groove Secures Lower U-Levers to Carriage Mounting Brackets	4		
Same as Part#	100-8601-001	.516ID x 1.25OD x .062Thk Nylon Washer Seats between Long Mounting Brackets & U-Lever/Caster Bases			
Same as Part#	100-8012-001	12mm External Snap Ring, Din 471 Between Lower U-Levers & Caster Carriage Brackets & Lock Nuts	4		
<b>CASTER &amp; WHEEL MOUNTING HARDWARE</b>					
Same as Part#	100-6744-001	7/16-14 Full Hex Head Nylon Lock Nut Secures Casters to Head & Foot Caster Bases/U-levers Crossbraces	4		
Same as Part#	100-5438-032	3/8-16 x 2.25" Hex Head Cap Screw Secures Rubber Wheels to Head & Foot Caster Bases/U-levers	4		
Same as Part#	100-6738-003	3/8-16 Thin Nylon Jam Lock Nut Secures Cap Screws for Wheels to Head & Foot Bases/U-levers	4		

CONTROL BOX AND MOTOR MOUNTING HARDWARE					
Purchase #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
HEAD MOTOR & CONTROL BOX MOUNTING HARDWARE					
-	Motor Clip	Motor Clip Comes with Control Box 999-0779-300 Secures Control Box to Head Motor Casing	1		
-	999-0775-129G	Head Actuator Bracket Attaches to Long Mounting Bracket with Carriage Bolts	1		
Same as Part#	999-0778-498	12mm Dia. x 98mm Clevis Pin with hole Secures Head Motor/Actuator to Head Actuator Bracket	1		
Same as Part#	999-0775-105	.504"ID x .62OD x .813Thk Nylon Spacer Seats between Head Motor Bracket and Long Mounting Bracket, Inserted onto 98mm Clevis Pin	1		
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Head Motor Shaft to Head/Back Deck Bracket	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Head Motor Clevis Pins to Brackets	2		
Same as Part#	100-5712-002	#12 x .50" Phillips Head Self Tapping Screw Secures Power Cord Ground Wire to Main Frame	1		
FOOT MOTOR MOUNTING HARDWARE					
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Foot Motor End to Main Frame Crossbrace Bracket	1		
Same as Part#	999-0775-445	12mm Dia. x 45mm Clevis Pin with hole Secures Foot Motor Shaft to Knee Deck Bracket	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Foot Motor Clevis Pins to Brackets	2		
HI/LO MOTOR MOUNTING HARDWARE					
Same as Part#	999-0775-445	12mm Dia. x 45mm Clevis Pin with hole Secures Hi/Lo Actuator Motors to Pan Seat Brackets - hd & ft hi/los	2		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Clevis Pin at Head Motor/Actuator Shaft to Head Caster Base/U-lever Mounting Bracket	1		
Same as Part#	100-2001-006	1-9/16 Hair Pin Clip Secures Clevis Pin to Head Hi/Lo Actuator Shaft & Mounting Bracket	1		

MISCELLANEOUS HARDWARE & INSTRUCTION LITERATURE					
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	100-8600-015	1.00" Dia. x .25"Thk Duopreen Cushion Disk (Head deck)	2		
Same as Part#	100-0510-001	Pendant Alligator Clip	1		
Same as Part#	999-0816-190	Service/Instruction Manual for Intro 705 Bed	1		



# STANDARD INTRO 705 BED SERVICE PART DIAGRAM

Model 705; NT7174 Shown

