



# ZENITH III BED ORDER FORM

Form# 999-0808-185  
REV C 02/10

## Model 708

To order NEW beds, please contact a GF Health Products Sales Representative at 1-800-347-5678.

ZENITH III BEDS AVAILABLE			
Order #	Description	Qty Ordered	Order Date
<b>ZG Series with Two Standard Caster Base Assemblies</b>			
ZG7074	76" Grid Deck — Full electric with 1 Head, 1 Foot, 2 Hi/Lo Motors, and 1 J-Box		
ZG7174	80" Grid Deck — Full electric with 1 Head, 1 Foot, 2 Hi/Lo Motors, and 1 J-Box		
ZG7014	76" Pan Deck — Full electric with 1 Head, 1 Foot, 2 Hi/Lo Motors, and 1 J-Box		
ZG7114	80" Pan Deck — Full electric with 1 Head, 1 Foot, 2 Hi/Lo Motors, and 1 J-Box		
<b>ZZ Series with One Step Locking Caster Base Assembly at foot end and Caster Base with Keepers at head end</b>			
ZZ7074	76" Grid Deck — Full electric with 4 Motors, 1 J-Box, and Z-Lok Locking System		
ZZ7174	80" Grid Deck — Full electric with 4 Motors, 1 J-Box, and Z-Lok Locking System		
ZZ7014	76" Pan Deck — Full electric with 4 Motors, 1 J-Box, and Z-Lok Locking System		
ZZ7114	80" Pan Deck — Full electric with 4 Motors, 1 J-Box, and Z-Lok Locking System		
<b>ZZX Fast Series with One Step Locking Caster Carriage System and TWO J-Boxes/Control Boxes</b>			
ZZX7177	80" Grid Deck — Full electric with 4 Motors, 2 J-Boxes, and Z-Lok Locking System		
ZZX7017	80" Pan Deck — Full electric with 4 Motors, 2 J-Boxes, and Z-Lok Locking System		

**NOTE: Footboard with embedded Staff/Nurse Control is optional and is not part of the standard Zenith III Bed.**

## TABLE OF CONTENTS:

For new or replacement bed service parts please refer to the following pages.

<b>Main Bed Section</b> (Decking, u-levers, radial arm tube, caster carriage, casters, etc.)	Page 2
<b>Boards With Staff Control Assemblies</b> (Hd/Ft Boards, Staff Control, Cables, Shroud, etc.)	Page 3
<b>Bed End Components</b> (Standard boards, embedded footboards, mounting kits, etc.)	Page 3
<b>Wallsavers and Mattress Retainers</b>	Page 3
<b>Electronic Components</b> (Cables, motors, control boxes, hand controllers/pendants, etc)	Page 4
<b>Assist Devices</b> (Half rails, foot mounting bracket, assist bars, trapeze adaptors, trapeze supports)	Page 5
<b>Accessories</b> (4" extension kits, wide kits with boards, I.V. mount, traction adaptor)	Page 5
<b>Wide Deck Kits</b> (without boards) & <b>4" Extension Pan Kit Components</b>	Page 6
<b>Bed Hardware</b>	Page 7
<b>Control Box and Motor Mounting Hardware</b>	Page 8
<b>Economy Lock Out &amp; Battery Backup Options</b>	Page 9
<b>Replacement Casters for Standard and Locking Caster Bases</b>	Page 9
<b>Miscellaneous Hardware and Literature</b>	Page 9
<b>Standard Zenith III Bed Diagram</b>	Page 10



# ZENITH III BED SERVICE PART ORDER FORM

Model 708

Form# 999-0808-185  
REV C 02/10

The following list of service parts for the Zenith III Bed have been provided for your convenience. Please keep this list in a secure location and use as a reference when ordering new or replacement parts. A general bed diagram has been provided on the last page.

**To order new or replacement SERVICE PARTS, please contact a Basic American Medical Products Customer Service Representative at 1-800-365-2338 to confirm you are ordering the right part.**

MAIN BED SECTIONS		Please see back diagram (page 10) for general part identification.			
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
DECKING					
Same as Part#	999-0808-900G	Main Frame Weldment	1		
Same as Part#	999-0779-920G	Standard Head/Back Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0779-930G	Standard Foot Grid Deck Weldment (80")	1		
	999-0775-936G	Optional Foot Grid Deck Weldment (76")	(1)		
Same as Part#	999-0779-940G	Standard Knee Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0796-982G	Optional Head/Back Pan Deck Weldment (76" & 80")	(1)		
Same as Part#	999-0796-983G	Optional Foot Pan Deck Weldment (80")	(1)		
	999-0796-986G	Optional Foot Pan Deck Weldment (76")	(1)		
Same as Part#	999-0796-984G	Optional Knee Pan Deck Weldment (76" & 80")	(1)		
Same as Part#	999-0748-001	One Edema Foot Ratchet Positioner	2		
Same as Part#	999-0775-910	Leveler Mount Assembly	2		
HEAD & FOOT U-LEVERS & DRAG LINK					
Same as Part#	999-0808-913G	Lower Head U-Lever Weldment	1		
Same as Part#	999-0808-914G	Upper Head U-Lever Weldment	1		
Same as Part#	999-0808-903G	Lower Foot U-Lever Weldment	1		
Same as Part#	999-0808-904G	Upper Foot U-Lever Weldment	1		
Same as Part#	999-0807-906G	Upper & Lower Drag Link Weldments	1 ea.		
Same as Part#	999-0808-012G	One Spanner Tube (Connect head and foot lower u-levers)	2		
CASTER BASE ASSEMBLIES					
Same as Part#	999-0808-980G	Optional 80" Pedal Lock Caster Base Assm. - Foot	1		
	999-0808-976G	Optional 76" Pedal Lock Caster Base Assm. - Foot	1		
	999-0808-958G	Caster Base with Keepers Assm. - Head end (76" & 80")	1		
Same as Part#	999-0808-955G	Standard Caster Base w/locking & non-locking casters	2		
CAPPING (if ordered separately)					
Same as Part#	999-0775-001	One Half Moon Cap (Upper & Lower U-Lever Ends)	10		
Same as Part#	999-0775-001	One Half Moon Cap (Four Each Caster Base)	8		
Same as Part#	100-4700-015	One 1.00" Sq. Tube End Cap (Main Frame Crossbraces)	4		
Same as Part#	100-4700-017	One 1.00" x 2.00" Rect. Tube End Cap (Drag Link Ends)	4		
Same as Part#	100-4700-021	One .75" x 1.50" 11GA. Rect. Tube End Cap (Ends of both Spanner Tubes)	4		

**HD/FT BOARDS WITH EMBEDDED STAFF CONTROL ASSEMBLIES**

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
ZL788053	-	Solar Oak Hd/Ft. Boards with Staff Control	1/1/1		
ZL788048	-	Med. Oak Hd/Ft. Boards with Staff Control	1/1/1		
ZL788080	-	Montana Walnut Hd/Ft. Boards with Staff Control	1/1/1		
ZL788029CCR	-	Windsor Mahogany Hd/Ft. Boards with Staff Control	1/1/1		
ZL788002CCR	-	Figured Mahogany Hd/Ft. Boards with Staff Control	1/1/1		
ZL788004CCR	-	Wild Cherry Hd/Ft. Boards with Staff Control	1/1/1		

**STANDARD BED END COMPONENTS**

Please see back diagram (page 10) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
<b>STANDARD HEAD- &amp; FOOTBOARDS WITHOUT CUTOUT:</b>					
Z53500	-	Solar Oak Hd/Ft Boards without staff control cutout	1/1		
Z48500	-	Med Oak Hd/Ft Boards without staff control cutout	1/1		
Z80500	-	Walnut Oak Hd/Ft Boards without staff control cutout	1/1		
Z29500CCR	-	Wind Mah CC Hd/Ft Boards without staff control cutout	1/1		
Z04500CCR	-	Wild Cher CC Hd/Ft Boards without staff control cutout	1/1		
Same as Part#	100-7025-004	1/4-20 UNC x .60 Long Threaded Inserts (pre-inserted into footboard and headboard)	8 (4 per bd)		
<b>STANDARD HEAD- &amp; FOOTBOARDS WITH STAFF CONTROL CUTOUT:</b>					
Z53650	-	Solar Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z48650	-	Med Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z80650	-	Walnut Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z29650CCR	-	Wind Mah CC Hd/Ft Boards cutout for emb. staff control	1/1		
Z04650CCR	-	Wild Cher CC Hd/Ft Boards cutout for emb. staff control	1/1		
Same as Part#	100-7025-004	1/4-20 UNC x .60 Long Threaded Inserts (pre-inserted into footboard and headboard)	8 (4 per bd)		
<b>STANDARD MOUNTING KIT 999-0775-000 INCLUDES:</b>					
Same as Part#	999-0746-035G	Angled Mounting Tubes (for mounting boards to frame)	4		
Same as Part#	100-4700-015	1" Square Tube End Caps (pre-inserted in mounting tubes)	4		
Same as Part#	100-5425-081	1/4-20 UNC x 1.50" Phillips Pan Head Machine Screws	8		
Same as Part#	100-7931-007	5/16 x 1.50" Clevis Pin with hole (mounting tubes to frame)	4		
Same as Part#	100-2001-006	1-9/16 Hair Pin Clip (mounting tubes to frame)	4		
Same as Part#	999-0775-195	Assembly Instruction Sheet	1		

**OTHER COMPONENTS**

Please see back diagram (page 10) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
ZZ6174	999-0775-180	Standard Wireform Head-end Wallsaver (Works with 76", 80", 84" & trapeze adaptor)	1		
ZA81700	999-0817-001	Optional Tubular Head-end Wallsaver Kit	1		
ZA81800	999-0818-001	Optional Tubular Side Wallsaver Kit	2		
Same as Part#	999-0814-000	Corner Wire Mattress Retainer Kit (set of two retainers)	2 in kit		
Same as Part#	999-0810-000	Corner Wire Mattress Retainers for Wide Kit Beds	2 in kit		
Same as Part#	999-0814-001	Single Wire Mattress Retainer	1		

**ELECTRICAL COMPONENTS**

Please see back diagram (page 10) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	999-0775-052	Head Actuator/Motor	1		
Same as Part#	999-0775-201	Head Actuator/Motor Cable (28") - Yellow O-Ring	1		
Same as Part#	999-0775-053	Foot/Knee Actuator/Motor	1		
Same as Part#	999-0775-202	Foot/Knee Actuator/Motor Cable (60") - Blue O-Ring	1		
Same as Part#	999-0778-051	One Standard Hi/Lo Actuator/Motor	2		
Same as Part#	999-0776-051	One Opt. Fast Hi/Lo Actuator/Motor for fast bed	2		
Same as Part#	999-0775-203	Long Lower Hi/Lo Motor Cable (60") - Red O-Ring	1		
Same as Part#	999-0808-204	Short Upper Hi/Lo Motor Cable (26") - Red O-Ring	1		
Same as Part#	999-0779-300	J-Box/Control Box	1		
Same as Part#	999-0775-208	J-Box/Control Box 3-Prong Power Cable with Ground	1		
Same as Part#	999-0776-205	Optional "T" Power Cable for Dual Control Boxes	1		
Same as Part#	999-0775-200	J-Box/Staff Control Ext. Cable (8") - Green O-Ring	1		
Same as Part#	999-0808-200	Standard Pendant "T" Cable - Green O-Ring For Beds with one control box - Has two male ends on "T"	1		
Same as Part#	999-0776-200	Optional Pendant "T" Cable - Green O-Ring For Fast Bed with two control boxes - Has two female ends on "T" instead of two male ends.	1		
Same as Part#	999-0808-299	End Cap for Standard Pendant "T" Cable	2		
Same as Part#	999-0808-305	Standard Pendant/Hand Controller	1		
Same as Part#	999-0791-000	Pendant Holster	1		
Same as Part#	999-0775-901	Optional Embedded Staff Control (with Cables & Bezel)	1		
Same as Part#	999-0775-006G	Wide Cable Shroud/Cover	1		
ZL788053	-	Staff Control Headboard/Footboard Assembly Kit (Includes Staff Control Switch Pad with Bezel and Cables, Shroud Instruction Sheet, Solar Oak Hd/Ft Boards, Attaching Hardware - all appropriate screws - and 4" Nylon Cable Ties)	1 kit		
ZL788000	-	Staff Control Assembly Kit (same as above - no boards)	1 Kit		
Same as Part#	100-1100-018	One 12" Nylon Cable Tie (to secure all cabling)	22		
	100-1100-050	One 17.5" Nylon Cable Tie (to secure control box)			
Standard Hi/Lo Actuator/Motor Kit (999-0808-950) with two control arm assemblies and hi/lo motors includes . . .					
Same as Part#	999-0778-051	Single Standard Hi/Lo Actuator/Motor	2		
Same as Part#	999-0808-004G	Single Central Arm Plate for Roller Arms	4		
Same as Part#	999-0779-103	Single Phenolic Wheel for Roller Arms	4		
Same as Part	100-8450-004	33/64"ID x 15/16"OD x 1/8"THK Washer	4		
Same as Part#	999-0808-487	12mm x 87mm Clevis Pin with Groove for Roller Arm	2		
Same as Part#	100-8012-001	12mm Snap Rings to secure Clevis Pins	2		
Same as Part#	999-0808-404	.50" x 1.50" Clevis Pin with Groove - Foot end	1		
Same as Part#	100-8050-002	Snap Ring - Foot-end Roller Assembly	1		
Optional Fast Hi/Lo Actuator/Motor Kit (999-0776-950) with two control arm assemblies and fast hi/lo motors includes . . .					
Same as Part#	999-0776-051	Single Opt. Fast Hi/Lo Actuator/Motor for fast bed	2		
Same as Part#	-	Everything else is the same as on Standard	-		

**ASSIST DEVICES**

NOTE: Rail "Kits" include both Quick Release Lanyard and Locking Bolt/Nut mounting hardware

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA78400 (both rails & hardware)	-	Left - Head Half Side Rails (76" & 80") Right - Head Half Side Rails (76" & 80")	1/1/4 for head end		
ZA78500 (includes both rails & two foot bracket kits)	-	Left - Foot Half Side Rails (80") with hardware Right - Foot Half Side Rails (80") with hardware Foot Bracket Kits with hardware	1/1/4 for 80" foot end		
ZA79400 (includes both rails & two foot bracket kits)	-	Left - Foot Half Side Rails (76") with hardware Right - Foot Half Side Rails (76") with hardware Foot Bracket Kits with hardware	1/1/4 for 76" foot end		
	999-0784-907G	Rt. Head (76"/80") or Lft. Foot (80") Half Rail only	1		
	999-0784-908G	Lft. Head (76"/80") or Rt. Foot (80") Half Rail only	1		
	999-0794-907G	Left Foot Half Rail for 76" only	1		
	999-0794-908G	Right Foot Half Rail for 76" only	1		
-	999-0786-995	Foot Bracket Hardware Bag with Bracket	1		
-	999-0775-025G	Foot Bracket only	1		
ZA80300 (includes hardware)		One Assist Bar with Mounting Hardware (Includes Quick Release Lanyard and Bolt/Nut Hardware)	1/2 for head end		
ZA80350 (includes 1 foot bracket)	-	One Assist Bar with mounting hardware, and Foot Bracket Kit with mounting hardware	1/2 for foot end		
-	999-0803-902G	One Assist Bar only (without hardware)	1		
-	999-0784-100	Quick Release Lanyard only	1		
ZA79000	-	Trapeze Support Adaptor Kit (Includes Trapeze Support and attaching hardware)	1 for head end		
ZA78100	-	Trapeze Bar Offset (must use with ZA79000)	1		

**ACCESSORIES**

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA71900	-	4" Standard Extension Foot Pan Kit	1 kit		
ZA82100		4" Extension Foot Pan Kit for Wide Deck	1 kit		
ZA76800	-	Traction Adaptor Kit	1		
ZA19000	554-0190-000	Brd. Mounted I.V. Sockets (includes bracket, cap & hdwr)	1		
WIDE KITS FOR ZENITH III BEDS					
ZA811000	-	42" Wide Kit for Rail Mount - No Boards			
WIDE KITS FOR ZENITH III BEDS WITH BOARDS & STAFF CONTROL ASSEMBLIES					
ZA811053	-	42" Wide Kit for Rail Mount - Solar Oak Bds w/Staff Control			
ZA811048	-	42" Wide Kit for Rail Mount - Med Oak Bds w/Staff Control			
ZA811080	-	42" Wide Kit for Rail Mount - Walnut Oak Bds w/Staff Control			
ZA811004CCR	-	42" Wide Kit for Rail Mount - CherCC Oak Bds w/Staff Control			
ZA811029CCR	-	42" Wide Kit for Rail Mount - WinCC Oak Bds w/Staff Control			

**42" WIDE DECK KIT**

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA811000	-	42" WIDE DECK KIT WITHOUT BOARDS	1		
Kit Includes the Following . . .					
-	999-0811-921G	Left-Hand Head Extension Weldment	1		
-	999-0811-922G	Right-Hand Head Extension Weldment			
-	999-0810-004G	Left-Hand Seat Pan Extension Weldment			
-	999-0810-005G	Right-Hand Seat Pan Extension Weldment			
-	999-0810-002G	Left-Hand Knee Extension Weldment			
-	999-0810-012G	Right-Hand Knee Extension Weldment			
-	999-0811-931G	Left-Hand Foot Extension Weldment			
-	999-0811-932G	Right-Hand Foot Extension Weldment			
-	999-0811-950	Hardware Box	1		
999-0811-950 Hardware Box Includes . . .					
-	999-0032-010	Short Self-Wrenching Nut (Head end pan assemblies to foot end pan assemblies)	2 1 each side		
-	100-5431-002	5/16-18 x .375" Hex Washer Head Machine Screw (Head end pan assemblies to foot end pan assemblies)	2 1 each side		
-	100-7931-008	5/16 x 3.50" Clevis Pin w/Hole (Secures pan sides to deck frame sides (4 RH and 4LH))	8 4 each side		
-	100-2001-006	1/9/16 Hair Pin Clip (Secures long clevis pins to deck frame)	8 4 each side		
-	999-0811-015	Wide Kit Instruction Sheet	1		

**EXTENSION PAN KITS**

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA71900	-	STANDARD 4" EXTENSION PAN KIT includes:	1		
-	999-0719-001G	Standard 4" Extension Pan - Grey	1		
-	999-0719-995	Hardware Bag	1		
Hardware Bag Includes . . .					
-	100-5231-001	5/16-18 x .75" Carriage Bolt (Attaches pan to foot deck from top)	2 1 each side		
-	100-6931-001	5/16-18 Wing Nut (Secures carriage bolts attaching pan to foot deck)	2 1 each side		
-	100-7931-006	5/16 x 1.75" Clevis Pin			
-	999-0719-115 (p1)	4" Expansion Pan Instruction Sheet	1		
Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA821000	-	WIDE DECK 4" EXTENSION PAN KIT	1		
Kit Includes the Following . . .					
-	999-0821-001G	Wide Deck 4" Extension Pan - Grey	1		
-	999-0821-995	Hardware Bag	1		
Hardware Bag Includes . . .					
-	100-5225-009	1/4-18 x .75" Carriage Bolt (Attaches pan to foot deck from top)	2 1 each side		
-	100-8525-001	1/4-18 Spring Lock Washer (Sits on carriage bolt between deck bottom and pan knob)	2 1 each side		
-	999-0810-009	Black Pan Knob (Secures carriage bolts attaching pan to foot deck)	2 1 each side		
-	999-0719-115 (p2)	Wide Deck 4" Expansion Pan Instruction Sheet	1		

**BED HARDWARE**

Purchase #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
<b>DECK TO FRAME MOUNTING HARDWARE</b>					
Same as Part#	100-5431-002	5/16-18 x .375" Hex Washer Head Machine Screw Secures Head & Knee Decks to Pan Seat and Knee to Foot Deck,	6		
Same as Part#	999-0032-010	Short Self-Wrenching Nut Secures Head & Knee Decks to Pan Seat and Knee to Foot Deck	6		
Same as Part#	100-7950-004	.50" x 2.00" Clevis Pin with hole Secures Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-7518-001	7/16-.50" Circle Cotter Secures Clevis Pin for Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-5431-023	5/16-18 x 2.55" Hex Head Machine Screw Secures Edema Ratchets to Frame Leveler Mounting Weldments	2		
Same as Part#	100-5231-005	5/16-18 x 1.75" Carriage Bolt Secures Leveler Strap Assembly to Main Frame Rail	2		
Same as Part#	100-6731-006	5/16-18 Thin Nylon Lock Nut Secures Carriage Bolts for Leveler Straps to Frame Leveler Mounting Weldment & Ratchet Shoulder Bolts to Foot Deck	6		
Same as Part#	999-0633-100	5/16-18 x .50" Hex Head Shoulder Bolt Secures Edema Ratchets to Foot Deck Mounting Weldments	2		
Same as Part#	100-8410-004	.437ID x .625OD x .093THK Spacer Between Foot Deck Edema Ratchets and Nylon Lock Nut	2		
<b>U-LEVERS TO MAIN FRAME AND DRAG LINK MOUNTING HARDWARE</b>					
Same as Part#	999-0779-401	.50" x 1.50" Clevis Pin with groove Secures Upper U-Levers to Main Frame Brackets	2		
Same as Part#	100-8050-002	.50" Dia. External Snap Ring Secures Upper U-Lever Clevis Pins to Main Frame Brackets	2		
Same as Part#	100-8601-001	.516"ID x 1.25"OD x .062"THK Nylon Washer Between Upper U-Levers to Main Frame Brackets	4		
Same as Part#	100-7950-005	.50" x 1.75" Clevis Pin with hole Secures Upper & Lower Drag Links to U-Levers	4		
Same as Part#	100-7950-006	.50" x 3.50" Clevis Pin with hole Secures Spanner Tubes to both Upper & Lower U-Levers	4		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures U-lever and Drag Link Clevis Pins	8		
Same as Part#	806-0277-008	Long Self-Wrenching Nut Secures Upper Head U-lever to Main Frame Brackets	2		
Same as Part#	100-5431-002	5/16-18 x .375 Hex Washer Head Machine Screw Secures Upper Head U-lever to Main Frame Brackets	2		
<b>CASTER CARRIAGE MOUNTING HARDWARE</b>					
Same as Part#	999-0779-401	.50" x 1.50" Clevis Pin with groove Secures Lower U-Levers to Caster Base Mounting Brackets	4		
Same as Part#	100-8050-002	.50" Dia. External Snap Ring Secures Lower U-Lever Clevis Pins to Caster Base Brackets	4		
Same as Part#	100-6744-001	7/16-14 Full Nylon Lock Nut - Standard Carriage Secures Caster Stems to Caster Carriage Mounting Tubes	8		
<b>CASTER KEEPER HARDWARE - LOCKING CARRIAGE ONLY</b>					
Same as Part#	999-0775-150G	Inner Keeper 2 used on Standard Carriage Caster (head end)	2		
Same as Part#	999-0775-151G	Outer Keeper 2 used on Standard Carriage Caster (head end)	2		
Same as Part#	100-5425-059	1/4-20 x 2.50" Hex Head Cap Screw Used to secure Inner Keeper to Outer Keeper	2		
Same as Part#	100-6725-003	1/4-20 Thin Nylon Lock Nut Secures Hex head Cap Screw to Inner Keeper and Outer Keeper	2		

## CONTROL BOX AND MOTOR MOUNTING HARDWARE

Purchase #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
CONTROL BOX MOUNTING HARDWARE					
-	Cannot be ordered separately	Control Box Clip (comes with control box) Secures Control Box to Head Motor Chassis	1 per box		
-	Cannot be ordered separately	Control Box Cable Clip (comes with control box) Secures Cables to Control Box ports	1 per box		
Same as Part#	999-0776-028G	Front Mounting Bracket for Control Box Secures to Frame Cross Brace - Control Box slides on	1 per box		
Same as Part#	999-0775-027G	Side Mounting Bracket for Control Box Secures to inside of Main Frame Rail using existing Bolts & Nylon Lock Nuts used to attach Leveler Strap Assembly to Main Frame. Control Box is secured to Bracket with 17.5" Cable Tie.	1 per box		
Same as Part#	100-5225-006	1/4-20 x 1.50" Bronze Carriage Bolt Secures Front Mounting Bracket to Main Frame Cross Brace	2 per bracket		
Same as Part#	100-6725-002	1/4-20 Full Height Nylon Lock Nut Secures Front Mounting Bracket to Main Frame Cross Brace	2 per bracket		
Same as Part#	100-1100-050	17.5" Nylon Cable Tie Secures Control Box to Head Motor Chassis	1 per box		
Same as Part#	100-5712-002	#12 x .50" Phillips Head Self-Tapping Screw Secures Ground Wire to Main Frame Rail	1 only		
HEAD MOTOR MOUNTING HARDWARE					
Same as Part#	999-0778-498	12mm Dia. x 98mm Clevis Pin with hole Secures Head Motor End to Main Frame Mounting Bracket	1		
Same as Part#	999-0775-105	.50ID x .62OD x .813THK. Spacer Inside Main Frame Mounting Bracket around Clevis Pin	1		
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Head Motor Shaft to Head Deck Bracket	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Head Deck & Motor End Mounting Clevis Pins	2		
FOOT MOTOR MOUNTING HARDWARE					
Same as Part#	999-0775-445	12mm Dia. x 45mm Clevis Pin with hole Secures Foot Motor End to Main Frame Crossbrace Bracket	1		
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Foot Motor Shaft to Knee Deck Bracket	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Foot Mounting Clevis Pins	2		
HI/LO MOTOR ASSEMBLY MOUNTING HARDWARE (For Standard and Fast Beds)					
Same as Part#	999-0775-445	12mm Dia. x 45mm Clevis Pin with hole Secures Hi/Lo Motor Ends to Upper & Lower Drag Link Brackets	2		
Same as Part#	100-2001-006	1-9/16 Hair Pin Clip Secures Hi/Lo Motor Ends/Upper & Lower Drag Link Clevis Pins	2		
Same as Part#	999-7950-004	.50" x 2.00" Clevis Pin with hole Secures Lower Hi/Lo Motor Control Arms to Head Lower Drag Link	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Lower Hi/Lo Control Arm Clevis Pins to Lower Drag Link	1		
Same as Part#	999-0808-404	.50" x 2.00" Clevis Pin with groove Secures Upper Hi/Lo Motor Control Arms to Foot Upper Drag Link	1		
Same as Part#	100-8050-002	.50" Dia. External Snap Ring Secures Upper Hi/Lo Control Arm Clevis Pins to Upper Drag Link	1		

**ECONOMY LOCK OUT AND BACKUP BATTERY CHARGER**

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
ZA78900	-	Economy Lockout Box (Includes mounting bracket & hardware)	1		
Same as Part#	999-0711-002	Battery Charger/Transformer (order with battery pack)	1		
Same as Part#	999-0711-003	Backup Battery Pack (order with charger/transformer)	1		

**REPLACEMENT CASTERS**

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	100-6330-061	Single Non-Locking 3" Twin Casters for Std. Caster Base	2 per base		
Same as Part#	100-6330-062	Single Locking 3" Twin Casters for Std. Caster Base	2 per base		
Same as Part#	100-6330-061	Single Non-Locking 3" Twin Casters for Z-Lok Caster bases	4 per base		

**MISCELLANEOUS HARDWARE & INSTRUCTION LITERATURE**

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	100-8600-015	1.00" Dia. x .25"Thk Duopreen Cushion Disk (head deck)	2		
Same as Part#	100-4550-010	Button Stops (bottom of main frame brackets)	4		
Same as Part#	100-0510-001	Pendant Alligator Clip (for pendant cable)	1		
Same as Part#	999-0808-190	Service/Instruction Manual for Zenith III Bed	1		



# STANDARD ZENITH III BED SERVICE PART DIAGRAM

## Model 708 Series

