


Five-Leg Deluxe I.V. Stand

- Chrome-plated, steel I.V. pole
- Can be easily converted from a 2-hook to a 4-hook pole
- Five-leg base offers maximum stability
- Swivel ball-bearing casters provide smooth rolling on all surfaces
- Base is molded of truss-reinforced plastic
- Sure-Grip clutch nut provides a full range of height positions from 47" to 81"
- Designed to meet the requirements of HCPCS code: E0776

GF7050-1 2-hook Top 1 ea

GF7051-1 4-hook Top 1 ea



GRAHAM-FIELD

www.grahamfield.com