

## Non-Rebreather Masks

- Used in medical emergencies for those with physical trauma that require high-flow oxygen, but do not require breathing assistance
- Designed to cover both the nose and mouth and attaches with the use of an elastic strap around the head for secure fit
- Comes with attached reservoir bag and 7' tubing
- Available with and without safety vent

<b>GF64159</b>	Adult without safety vent	50/cs
<b>GF64160</b>	Adult with safety vent	50/cs
<b>GF64190</b>	Pediatric with safety vent	50/cs



**GF64159**