



Basic Wedge Cushion

Three layers provide pressure distribution and weight equalization:

- 1. High density foam base
- 2. Non-aqueous specialized visco-elastic gel core layer reduces tissue interface pressure at critical points
 - Proprietary gel is a flexible solid that will not leak
 - Gel is placed at three pressure zones
- 3. Convoluted foam topper

Cover features

- Fluid-proof stretch top cover, reduces shear/friction
- Antimicrobial
- · Flame retardant
- · Removable and washable
- · Non-skid base

4" in front tapers to 3 ¼" in rear Strap to secure cushion firmly to the chair 350 lb maximum weight capacity 1-year limited warranty

Designed to meet the requirements of HCPCS Code: E2601

 8220166
 16" x 16" x 3 ¼ - 4"

 8220186
 18" x 16" x 3 ¼ - 4"

 8220188
 18" x 18" x 3 ¼ - 4"

 8220206
 20" x 16" x 3 ¼ - 4"

 8220208
 20" x 18" x 3 ¼ - 4"

Designed to meet the requirements of HCPCS Code: E2602

8220228 22" x 18" x 3 ¼ - 4" **8220248** 24" x 18" x 3 ¼ - 4"









