

Tenderflo® Water Overlay

- Five individual baffles allow even movement of fluid to provide the right balance of pressure reduction and support
- Two sealed seams provide maximum strength and durability build-up
- Heavy duty vinyl cleans easily with mild detergent
- Fill lines and foot placement lines can be used as guides for patient positioning
- Thermometer is included to monitor temperature
- Holds 15-20 gallons
- 225 lb maximum weight capacity
- Designed to meet the requirements of HCPCS code: E0198
- 90-day limited warranty

AQTF-2035 68" x 32" x 3 1/4" 1 ea



GRAHAM-FIELD

www.grahamfield.com

© 2011 GF Health Products, Inc.
Atlanta, Georgia 30360 GF1100053RevA11