

Gel-Lite III Gel Overlay

Gel-Lite III Gel Overlay

- A uniform 2" layer of gel covers 80% of the surface area of overlay and provides the right balance of pressure relief and support
- Features a lightweight gel formula that reduces overall product weight
- The individual cell construction of the gel prevents pooling and gel migration so that the patient will not bottom out
- Fluid-resistant polyurethane-coated nylon cover is durable, antibacterial, antifungal, and easily cleaned
- Four elastic straps secure the corners snugly to the mattress
- Conforms to any adjustable bed
- 300 lb maximum weight capacity
- Designed to meet the requirements of HCPCS code: E0185
- 1 year limited warranty

HDC-0072Z 78" x 34" x 3" 1 ea





www.grahamfield.com

© 2012 GF Health Products, Inc. Atlanta, Georgia 30360 GF1200059RevA12