




TRAVELER[®]

PARTS CATALOG

(Standard Model Only)



CONTENTS

FRAME—FIXED ARM, CROSSBRACE AND ARM PAD.....	3
FRAME—DETACHABLE ARM AND CROSSBRACE	4
ARM	5
REAR WHEEL	6
FOOTREST	7
ELEVATING LEGREST	7
UPHOLSTERY	8
CASTER.....	8
WHEEL LOCK	9
ACCESSORIES	9

TO ORDER

Tel: 770-368-4700 • Fax: 770-368-2386


Returns: Parts may not be returned without prior authorization from GF Health Products, Inc. Custom parts cannot be returned for credit. Please refer to GF Health Products, Inc. [Terms and Conditions](#) for further information.

GF Health Products, Inc. is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this catalog are based on the latest product information available at the time of printing. The most current product information, including the most current revision of this catalog, can be found online at www.grahamfield.com.


Everest & Jennings®, E&J®, Traveler®, and Graham-Field® are registered trademarks of GF Health Products, Inc.

© 2012, GF Health Products, Inc., All Rights Reserved

FRAME—FIXED ARM, CROSSBRACE AND ARM PAD




- 1 Sideframe, fixed arm; includes side panel and caster stem bearings..... 9P188217**
- 2 Crossbrace, pair; includes center bolt and nut**
 16" W 90761800
 18" W 90761801
 20" W 90761803
- 3 Inner seat rail, pair 9P516410**
- 4 Bottom inner rail kit; includes bottom inner rail, 2 each; screw and nut that hold bottom inner rail, 4 each; and sleeve, 2 each..... 9P2731A1**
- 5 Grip, black, 10 each..... 9BD82845**
- 6 Seat rail guide, 4 each..... 9PE25300**
- 7 Tipping sleeve, 10 each 9B279690**




Arm pad, pair; includes mounting screws

- Deep blue..... 9P190547
 Jet black 9P190539
 Vintage burgundy 9P1905P5

FRAME—DETACHABLE ARM AND CROSSBRACE



- | | | |
|----------|---|-----------------|
| 1 | Sideframe, detachable arm; includes | |
| | caster stem bearings | 9P250015 |
| 2 | Crossbrace, pair; includes center bolt and nut | |
| | 16"W..... | 90761800 |
| | 18"W..... | 90761801 |
| | 20"W..... | 90761803 |
| 3 | Bottom inner rail kit; includes | |
| | bottom inner rail, 2 each; screw | |
| | and nut that hold bottom inner rail, | |
| | 4 each; and sleeve, 2 each | 9P2731A1 |
| 4 | Inner seat rail, pair; includes plug..... | 9P516321 |
| 5 | Grip, black, 10 each | 9BD82845 |
| 6 | Seat rail guide 10 each | 9BE25320 |
| 7 | Tipping sleeve 10 each | 9B279690 |
| 8 | Seat slide kit, pair | |
| | 16" and 18" W..... | 9P06A002 |
| | 20" W..... | 9P064100 |




1 Arm, detachable fixed height, pair.
Includes side panel and spring lock;
does not include arm pad

Desk length 9P258353
Full length 9P232353

2 Spring lock, 10 each 9B232870

Arm pad: See below




1 Arm, detachable adjustable height, pair.
Includes side panel and spring locks;
does not include arm pad

Desk length 920000A5
Full length 920000A6

2 Spring lock, upper, 10 each 9BB32890

3 Spring lock, lower, 10 each 9B232870


Arm pad: See below



Arm pad, pair; includes mounting screws

Color	Desk length	Full length
Deep blue.....	9P190947	9P190547
Jet black.....	9P190939	9P190539
Vintage burgundy	9P1909P5.....	9P1905P5

REAR WHEEL



1 Rear wheel: size, tire, handrim options

24" x 1" black composite wheel, solid rubber tire, integral handrim	90-011B1
24" x 1" black composite wheel, solid rubber tire, chrome handrim	90-01J01
24" x 1 3/8" black composite wheel, pneumatic tire, chrome handrim	9000B008
24" x 1 3/8" black composite wheel, airless tire, chrome handrim	9000B016

Replacement parts included in item 1:

2 Hubcap, black, 10 each


1" tire	9B045035
1 3/8" tire	9B031366

3 Bearing, 10 each

1" tire	9B007700
1 3/8" tire	
Current wheel *	9B007710
Previous wheel **	9B850045

* without I.D. adapter bushing

** with I.D. adapter bushing



Replacement parts not included in item 1:

4 Axle and nut, 10 each

Fixed arm

1" tire	9B-07765
1 3/8" tire	9B-47320


Detachable arm

1" tire	9B-07845
1 3/8" tire	9B-07770

5 Spacer, 10 each.....


Integral handrim (3/16")	9B007360
1 3/8" tire (1/4")	9B007280

FOOTREST



Footrest, detachable and swingaway, pair (left shown). Includes right and left hanger with release and #1 aluminum footplate with extension.90761750

ELEVATING LEGREST



Elevating legrest, complete, pair. Detachable and swingaway. Includes right and left hanger with release, deep blue legrest panel, and #1 aluminum footplate with extension.

Low-pivot (shown)90761751
High-pivot90761755

1 Elevating legrest, pair. Detachable and swingaway. Includes right and left hanger with release and #1 aluminum footplate with extension. Does not include legrest panel.


Low-pivot90761758
High-pivot90761759

2 Legrest panel, pair; includes mounting hardware

Deep blue.....9P467247
Jet black.....9P467239
Vintage burgundy9P4672P5

3 #1 footplate with extension, pair.....9P426J03

UPHOLSTERY




Back upholstery (includes mounting hardware)

	16"W	18"W	20"W
Deep blue	99-2-247	99-4-047	99-4-247
Jet black	99-2-239	99-4-039	99-4-239
Vintage burgundy	99-2-2P5	99-4-0P5	99-4-2P5

Seat upholstery (includes mounting hardware)

	16"W	18"W	20"W
Deep blue	99-1-247	99-3-047	99-3-247
Jet black	99-1-239	99-3-039	99-3-239
Vintage burgundy	99-1-2P5	99-3-0P5	99-3-2P5

CASTER



Caster, complete, 8" x 1" solid tire;
includes caster with axle bolt,
axle nut and fork..... 90051910

1 Caster stem bearing, 10 each..... 9B291976

2 Caster fork, steel, with stem nut..... 92050355

3 Dust cap, black, 10 each..... 9B291975


4 Caster wheel, black molded, each

8" x 1" solid 90010015

8" x 1 3/4" pneumatic..... 90-51625

8" x 1 3/4" semi-pneumatic 90-51605

WHEEL LOCK



Wheel lock, pair; includes rubber tip

Push to lock

Fixed arms 9P761063
Detachable arms 9P761061

Pull to lock

Fixed arms 9P761067
Detachable arms 9P761065

Rubber tip, 10 each 9BD82746

Wheel lock extension, 10 1/2", pair 9A000026

ACCESSORIES

Accessories (not shown)

1 Positioning Strap

Auto-style push button buckle ... 92133040
Hook and loop style..... 90133115

2 O2-IV holder combo..... B10320E

3 IV holder, telescopic 90361110

4 O2 cylinder carrier GFM905-1

5 IV holder, one-hook..... GF1844W

6 Anti-tipper, pair..... 9A000020

7 Legrest anti-removal cap, 10 each .. 9B431140

8 Heel loop, pair 9A000002

9 Anti-fold bar..... 92523850



GF Health Products, Inc. ("Graham-Field")
2935 Northeast Parkway
Atlanta, Georgia 30360 USA

tel: 770-368-4700

fax: 770-368-2386



GRAHAM-FIELD

www.grahamfield.com

Traveler Parts Catalog
GF0700120RevC12