

## Rehab Shower Commode

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or used as a portable, self-contained commode.

- Completely removable full length arms for ease of transfer
- Swingaway, detachable footrests permit close approaches to beds, toilets and wheelchairs
- Recessed seat and frame side contours allow easy access for toileting activities while seated
- Available with durable 24" one piece composite wheels with integral handrims and solid rubber tires or 5" rear casters
- Optional plastic coated handrims for improved grip
- Removable, hygienic plastic pail and lid are included
- Cushioned vinyl seat is sealed against moisture penetration and vented for optimal foam compression
- Non-corrosive composite footplates
- Maximum weight capacity 250 lb - EVENLY DISTRIBUTED
- 2 year limited warranty on sideframes and crossbraces
- Chair comes standard with oval opening seat. C-shaped and solid seat are optional items



Optional C-shaped and solid seat available



Rehab Shower Commode		Dimensions and Weights							
		Seat			Back	Overall			Wt. *
12022010	14 1/2" LoBack	17"	18"	23 1/2"	14 1/2"	40"	38"	25"	43 lb
12022011	14 1/2" LoBack, Caster Pin Locks	17"	18"	23 1/2"	14 1/2"	40"	38"	25"	43 lb
12022013	14 1/2" LoBack with 5" Rear Locking Caster	17"	18"	23 1/2"	14 1/2"	40"	38"	25"	43 lb
12022110	17 1/2" HiBack	17"	18"	23 1/2"	17 1/2"	40"	41"	25"	45 lb
12022111	17 1/2" HiBack, Caster Pin Locks	17"	18"	23 1/2"	17 1/2"	40"	41"	25"	45 lb

All dimensions are ±1/4". \*Without Front Rigging.

