

NaviKnee Walker

The Lumex Knee Walker is ideal for individuals recovering from sprains, foot surgery, amputation or any foot or ankle injury where the user cannot bear weight on the injured leg. It can also be a more comfortable alternative to crutches where extended use can lead to underarm and/or back pain. The Lumex NaviKnee Walker can assist with a speedy recovery and help an individual maintain a more independent and active lifestyle.

Features:

- Lightweight, durable aluminum frame weighs less than 18 lbs
- Turning system allows for easy maneuverability & turning radius
- Deluxe dual braking system allows for quick and extended stopping
- Knee pad is adjustable both in height and distance from handle bars for a comfortable fit dependant on user's arm length
- Adjustable handle secured with easy turn knob to secure accommodates various user heights
- Deluxe polyurethane knee pad is designed for ultimate comfort
- Large, 8" casters with precision ball bearings are ideal for indoor/outdoor use
- Easy folding mechanism makes transport and storage simple
- Lifetime (3 Years) Limited Warranty
- Designed to meet the requirements of E0118

LX3000	Silver	1 ea
LX3000-POUCH	Black	1 ea



Specifications

Overall Width	13"
Overall Length	30"
Minimum Handle Height	34"
Maximum Handle Height	39"
Approx Patient Height.....	5' 1" – 6' 2"
Knee Pad Dimensions	14" x 8" x 4" (LxWxH)
Knee Pad Minimum Height.....	21"
Knee Pad Maximum Height.....	24"
Maximum Weight Capacity	300 lb EVENLY DISTRIBUTED
Product Weight	17.75 lb (8.05 kg)



GRAHAM-FIELD

www.grahamfield.com

© 2013 GF Health Products, Inc.
Atlanta, Georgia 30360 GF1200127RevA13