

## Digital Blood Pressure Monitor with Manual Inflation

**Latex Free**

- Reliable readings with simple operation
- Perfect for home use
- Manual inflation/auto deflation
- Offers easy-to-read, large digital display with a self-adjusting cuff for one-handed use
- Accurately measures systolic/diastolic pressures and pulse
- Automatic power-off function
- Compact size
- Latex-Free
- Instruction manual included
- One year warranty

**707A** Adult 1 ea

**707A C** Small Adult 1 ea

**707A X** Large Adult 1 ea

Cuff Class	Limb Range (Inches)	Cuff Dimensions (Inches)*
Newborn	3.75" – 5.25"	7" W x 2" H
Infant	5" – 7.5"	10.5" W x 3" H
Child	7" – 10.5"	13.5" W x 4" H
Adult	9.5" – 15"	20.5" W x 5.5" H
LG Adult	13" – 19"	25" W x 6.5" H
Thigh	17" – 26"	29.5" W x 8.5" H

\* These are subjective measurements and can be vary  $\pm 0.25"$

