


Silver Nitrate Applicators

- Applications: skin or mucous membrane cauterization, removing warts and granulated tissue
- 75% silver nitrate and 25% potassium nitrate
- Silver nitrate tips are soluble in all secretions
- Control strength by dipping in sterile, distilled water
- Rigid wooden applicator for precision placement
- Convenient tube vial
- Made in USA, GF Health Products, Inc. Rhode Island Facility
- ISO 13485:2003 Certified Facility

U.S. Federal Law prohibits dispensing without prescription.

Distributors (USA) must have an active state wholesale drug license to purchase Grafco® Silver Nitrate Applicators. Contact our Regulatory Affairs Dept to apply for Authorized Distributor of Record (ADR) credentials at 1800-556-6508.

1590	6", Double Dipped	NDC 12165-100-03	50mg active ingredients/applicator	100/vial, 12vials/box
1592	12", Double Dipped	NDC 12165-200-03	70mg active ingredients/applicator	100/vial
1593	18", Double Dipped	NDC 12165-300-01	98mg active ingredients/applicator	1 ea



1590



1592



1593