




***Traveler*[®] L4**

PARTS CATALOG



CONTENTS

SIDEFAME AND CROSSBRACE 3

ARM: FLIP-BACK FIXED-HEIGHT 4

ARM: FLIP-BACK ADJUSTABLE-HEIGHT 4

REAR WHEEL 5

CASTER..... 5

FRONT RIGGING, DETACHABLE 6

UPHOLSTERY 7

ACCESSORIES 7

TO ORDER

Tel: 770-368-4700 • Fax: 770-368-2386


Returns: Parts may not be returned without prior authorization from GF Health Products, Inc. Custom parts cannot be returned for credit. Please refer to GF Health Products, Inc. Terms and Conditions for further information.

GF Health Products, Inc. is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this catalog are based on the latest product information available at the time of printing. The most current product information, including the most current revision of this catalog, can be found online at www.grahamfield.com.

Everest & Jennings, E&J, Traveler, and Graham-Field are registered trademarks of GF Health Products, Inc.

© 2009, GF Health Products, Inc., All Rights Reserved

SIDEFRAME AND CROSSBRACE



1 Sideframe, all widths

16" deep Right.....	90763201
16" deep Left.....	90763202
18" deep Right.....	90763203
18" deep Left.....	90763204

2 Crossbrace, pair (Includes center bolt and nut)

16" W x 16" deep.....	907632C1
18" W x 16" deep.....	907632C2
20" W x 16" deep.....	907632C3
18" W x 18" deep.....	907632E2
20" W x 18" deep.....	907632E3

3 Crossbrace link, set of 4 (includes mounting hardware)

16" W x 16" deep.....	907632D1
18" W x 16" deep.....	907632D1
20" W x 16" deep.....	907632D3
18" W x 18" deep.....	907632D1
20" W x 18" deep.....	907632D3

4 Crossbrace saddle bushing and gasket, pair (Includes mounting hardware)

90763205

5 Grip (10 each)

90763005

6 Seat rail cradle, pair (includes 2 each screw, washer and nut)

90763206

7 Tipping sleeve (10 each)

90763006

8 Backpost, knob and screw, pair

90763207

9 Seat rail and arm cradle, pair (includes 2 each screw, washer and nut)

90763208


10 Mounting ring (10 each)

90763209

11 Arm bracket, pair (includes 2 each screw, washer and nut)

90763210

ARM: FLIP-BACK FIXED-HEIGHT



- 1** Arm, flip-back fixed-height, pair (includes right and left arm with mounting hardware, side panel, and latch; does not include arm pad)

Desk length (shown)... 90763220

Full length..... 90763221


- 2** Arm pad, pair (includes two black arm pads and mounting screws)

Desk length..... 90763022

Full length..... 90763023

- 3** Arm latch, pair, with hardware 90763223

ARM: FLIP-BACK ADJUSTABLE-HEIGHT



- 1** Arm, flip-back adjustable-height, pair (includes right and left arm with mounting hardware, side panel, and latch; does not include arm pad)

Desk length (shown)... 90763225

Full length..... 90763226


- 2** Arm pad, pair (includes two black arm pads and mounting screws)

Desk length..... 90763022

Full length..... 90763023


- 3** Arm latch, pair, with hardware 90763223

REAR WHEEL




- | | |
|---|----------|
| 1 Rear wheel with handrim, 24" x 1", polyurethane tire (includes hub cap and bearing; does not include axle or nut) | 90763150 |
| 2 Hub cap (10 each)..... | 90763051 |
| 3 Wheel bearing (10 each) | 90763052 |
| 4 Rear axle and nut (10 each)..... | 90763053 |
| 5 Quick-Release axle: for use with adjustable-height-arm wheelchair model <u>only</u> (10 each) | 90763252 |

CASTER




- | | |
|---|----------|
| Complete caster (fork and wheel), 8" x 1" solid rubber tire | 90763040 |
| 1 Caster wheel, axle and nut | 90763041 |
| 2 Caster fork, stem nut and dust cap | 90763042 |
| 3 Caster stem bearing (10 each) | 90763043 |
| 4 Caster dust cap (10 each)..... | 90763044 |
| 5 Caster axle and nut (10 each)..... | 90763045 |

FRONT RIGGING, DETACHABLE




Footrest, pair, all widths
Includes right and left hanger with release
and composite footplate with extension
and heel loop 90763030

Composite footplate, pair
(Fits all front rigging) 90763033




Low-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
composite footplate with extension, and
legrest panel (does not include
heel loop) 90763031

Elevating legrest panel kit
(Fits all elevating legrests)
Includes 2 black padded panels
and hardware 90763032



High-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
aluminum footplate with extension, and
legrest panel (does not include
heel loop) 90763435

UPHOLSTERY



Back upholstery, black nylon, and mounting hardware (includes replacement ties)

16" W	907632A1
18" W	907632A2
20" W	907632A3

Seat upholstery, black nylon, and mounting hardware

16" W x 16" D.....	907632B1
18" W x 16" D.....	907632B2
20" W x 16" D.....	907632B3
18" W x 18" D.....	907632F2
20" W x 18" D.....	907632F3

ACCESSORIES

Traveler L4 accessories (not shown)

1 Heel loop, pair.....	90763064
2 Wheel lock, pair (includes rubber tips)	
Push to lock.....	90763070
Pull to lock.....	90763071
3 Rubber tip, wheel lock (10 each).....	90763063
4 Wheel lock extension, pair	90763062
5 Positioning belt (metal buckle)	92133040
6 Anti-Fold bar	92523850
7 O2-IV holder	B10320F
8 O2 cylinder holder.....	B10325
9 Anti-Tippers, pair.....	90763230

Everest & Jennings®

GF Health Products, Inc. ("Graham-Field")
2935 Northeast Parkway
Atlanta, Georgia 30360 USA

tel: 770-368-4700

fax: 770-368-2386



GRAHAM-FIELD

www.grahamfield.com

Traveler L4 Parts Catalog
GF0900085RevG14