


Zenith Bed Series Accessories



ZA84200





Optional Half Length Assist Device and Assist Bars provide a sturdy and secure hand hold to assist residents in and out of bed. Assist devices are removable without tools and meet FDA Entrapment Guidelines.

- Move out of the way for lateral transfers.
- Pivot away from egress area.
- Provide toe clearance even in lowest position to prevent staff injury.
- Hand pendant mounts facing inward or outward and away from hand grab area.



Pivot Assist Bar ZA85400 (tool-less)
ZA87900 (bolt-on)



Fixed Assist Bar ZA85000 (tool-less)
ZA87800 (bolt-on)



For APS9000



For Zenith9000,
7000, and 5000

Pendant Hand Control and Holster


Easy operation contoured pendant is shock-resistant, splash-proof and durable. Holster conveniently secures to bedrails and assist devices for easy access.



Wide Deck Kit (for non-Slide-W-I-D-E beds)
Optional wide kit adjusts the bed sleep deck to 39" or 42" width for greater resident comfort. (Head and foot boards and staff control electronics not included.)

ZA858000 39" wide deck kit

ZA851000 42" wide deck kit



May be attached with bolts or quick release pins.



Trapeze Adapter

Optional trapeze adapter allows for easy attachment of trapeze system to headboard and within easy reach for residents to reposition themselves.

ZA79000 (for use with
ZA78100 Trapeze Bar)



Length Extension Kit

Transform your bed to 84" with this optional 4-inch sleep deck extension.

ZA85350 Standard 35" width
ZA85300 for 39" or 42" widths
ZA88300 for Slide-W-I-D-E beds



Battery Backup Kit for Zenith9000, 7000, and 5000

Provides the ability to perform bed functions in the event of power outage.

ZA84400 Battery Backup and Charger



Battery Backup for APS9000

Provides the ability to perform bed functions in the event of power outage.

999-0711-003 Battery backup

999-0711-002 Charger

Side Wall Saver Kit

(Set of 2)

Provides wall protection down the length of the bed
(For Zenith9000APS,
Zenith9000, and Zenith7000)

ZA86200



Plastic Pan Kit (field-installed)

Smooth surface for easy cleaning.
Can be removed.


ZA87500



Board Mounted IV Socket

(IV pole not
included)

ZA19000



Traction Adapter (set of 4)

ZA76800

Not shown