

Traveler® L4

Everest & Jennings introduces the latest in a long tradition of quality wheelchairs. The Traveler L4 is the newest high-strength lightweight wheelchair addition to the Traveler line, a brand respected for performance and durability.

The Traveler L4 is a high-strength lightweight K0004 Medicare coded wheelchair. The Traveler L4 has dual axles and multi-position caster forks for an adjustable seat height of 18 or 20 inches. The wheelchair comes standard with non-marking polyurethane tires, anti-tippers and a chart pocket.

The Traveler L4 is manufactured with the same commitment to quality Everest & Jennings is known for.

Standard Features

- Weighs less than 34 lbs. without front rigging
- Black nylon upholstery
- Chart pocket on back
- "Tool Free" adjustable front rigging
- Anti-tippers come standard
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Flip back arms
- Height adjustable back
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb (300 lb. on models with 18" seat depth)
- EVENLY DISTRIBUTED
- Meets Medicare code K0004

Warranty

Limited lifetime warranty on sideframe and crossbraces



Height Adjustable Arms
(on select models)



Composite Footplates




Chart Pocket



Flip Back Arms



Dual Rear Axle



Multi-position Caster



Adjustable Front Rigging



Industry Compatible
Front Rigging












Wheel Locks



GRAHAM-FIELD

www.grahamfield.com

© 2014 GF Health Products, Inc.
Atlanta, Georgia 30360 GF0900098RevD14

Traveler® L4		Dimensions and Weights								
		Seat			Back	Overall				Wt.*
										
3F020120	Flip Back Desk Arm, Swingaway Footrest	18"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	31.8 lbs
3F020130	Flip Back Desk Arm, Elevating Legrest	18"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	31.8 lbs
3F020140	Flip Back Full Arm, Swingaway Footrest	18"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	31.8 lbs
3F020150	Flip Back Full Arm, Elevating Legrest	18"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	31.8 lbs
3F020160	Flip Back Height Adjustable Desk Arms, Swingaway Footrest, Quick-Release Wheels	18"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	31.8 lbs
3F020170	Flip Back Height Adjustable Desk Arms, Elevating Legrest, Quick-Release Wheels	18"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	31.8 lbs
3F020220	Flip Back Desk Arm, Swingaway Footrest	16"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	23"	11.5"	31.6 lbs
3F020230	Flip Back Desk Arm, Elevating Legrest	16"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	23"	11.5"	31.6 lbs
3F020250	Flip Back Full Arm, Elevating Legrest	16"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	23"	11.5"	31.6 lbs
3F020260	Flip Back Height Adjustable Desk Arms, Swingaway Footrest, Quick-Release Wheels	16"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	23"	11.5"	31.6 lbs
3F020270	Flip Back Height Adjustable Desk Arms, Elevating Legrest, Quick-Release Wheels	16"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	23"	11.5"	31.6 lbs
3F020320	Flip Back Desk Arm, Swingaway Footrest	20"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	27"	11.5"	32 lbs
3F020330	Flip Back Desk Arm, Elevating Legrest	20"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	27"	11.5"	32 lbs
3F020340	Flip Back Full Arm, Swingaway Footrest	20"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	27"	11.5"	32 lbs
3F020350	Flip Back Full Arm, Elevating Legrest	20"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	27"	11.5"	32 lbs
3F020360	Flip Back Height Adjustable Desk Arms, Swingaway Footrest, Quick-Release Wheels	20"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	27"	11.5"	32 lbs
3F020370	Flip Back Height Adjustable Desk Arms, Elevating Legrest, Quick-Release Wheels	20"	16"	18" or 20"	16.5" or 18.5"	31.5"	36"	27"	11.5"	32 lbs
3F030120	Flip Back Desk Arms, Swingaway Footrest	18"	18"	18" or 20"	16.5" or 18.5"	33"	36"	25"	11.5"	32 lbs
3F030130	Flip Back Desk Arms, Elevating Legrest	18"	18"	18" or 20"	16.5" or 18.5"	33"	36"	25"	11.5"	32 lbs
3F030140	Flip Back Full Arms, Swingaway Footrest	18"	18"	18" or 20"	16.5" or 18.5"	33"	36"	25"	11.5"	32 lbs
3F030150	Flip Back Full Arms, Elevating Legrest	18"	18"	18" or 20"	16.5" or 18.5"	33"	36"	25"	11.5"	32 lbs
3F030160	Flip Back Height Adjustable Desk Arms, Swingaway Footrest, Quick-Release Wheels	18"	18"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	32 lbs
3F030170	Flip Back Height Adjustable Desk Arms, Elevating Legrest, Quick-Release Wheels	18"	18"	18" or 20"	16.5" or 18.5"	31.5"	36"	25"	11.5"	32 lbs
3F030320	Flip Back Desk Arms, Swingaway Footrest	20"	18"	18" or 20"	16.5" or 18.5"	33"	36"	27"	11.5"	32.2 lbs
3F030330	Flip Back Desk Arms, Elevating Legrest	20"	18"	18" or 20"	16.5" or 18.5"	33"	36"	27"	11.5"	32.2 lbs
3F030340	Flip Back Full Arms, Swingaway Footrest	20"	18"	18" or 20"	16.5" or 18.5"	33"	36"	27"	11.5"	32.2 lbs
3F030350	Flip Back Full Arms, Elevating Legrest	20"	18"	18" or 20"	16.5" or 18.5"	33"	36"	27"	11.5"	32.2 lbs
3F030360	Flip Back Height Adjustable Desk Arms, Swingaway Footrest, Quick-Release Wheels	20"	18"	18" or 20"	16.5" or 18.5"	33"	36"	27"	11.5"	32.2 lbs
3F030370	Flip Back Height Adjustable Desk Arms, Elevating Legrest, Quick-Release Wheels	20"	18"	18" or 20"	16.5" or 18.5"	33"	36"	27"	11.5"	32.2 lbs

*Without Front Rigging
All dimensions are ±1/4"