

The Traveler LX incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction and industry standard front rigging. The premium vinyl upholstery is reinforced with an extra heavy duty inner liner. This wheelchair comes in a chrome finish and provides the most effective combination K0001/K0002 wheelchair available.



Standard Features

- Blue premium vinyl upholstery
- Chart pocket on back
- "Tool Free" adjustable front rigging
- Chrome finish provides a long lasting and easy to clean finish
- Industry compatible front rigging
- Gray non-marking tires
- Two position axle and caster fork provide 2" seat height and adjustment
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002

Optional Features

- Rear Mount Anti-tippers

Warranty

Two (2) year limited warranty



Traveler® LX		Dimensions and Weights								
		Seat			Back	Overall				Wt.*
3E020100	Fixed Full Arm, Swingaway Footrest	18"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	26 ½"	12"	44 lb
3E020110	Fixed Full Arm, Elevating Legrest	18"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	26 ½"	12"	44 lb
3E020120	Detachable Desk Arm, Swingaway Footrest	18"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	26 ½"	12"	44 lb
3E020130	Detachable Desk Arm, Elevating Legrest	18"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	26 ½"	12"	44 lb
3E020140	Detachable Full Arm, Swingaway Footrest	18"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	26 ½"	12"	44 lb
3E020220	Detachable Desk Arm, Swingaway Footrest	16"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	24 ½"	12"	44 lb
3E020320	Detachable Desk Arm, Swingaway Footrest	20"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	28 ½"	12"	44 lb
3E020330	Detachable Desk Arm, Elevating Legrest	20"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	28 ½"	12"	44 lb
3E020340	Detachable Full Arm, Swingaway Footrest	20"	16"	18" or 20"	15 ½"	31 ½"	35 ½"	28 ½"	12"	44 lb

*This product is not made with natural rubber latex

