

The Advantage Wheelchair provides comfort and durability at an affordable price. Designed as a low maintenance wheelchair, it offers an attractive silver vein durable hammertone steel finish.



Chart Pocket



Industry Compatible Front Rigging

Standard Features

- Available in 16", 18" and 20" seat widths
- Black padded nylon and vinyl upholstery available with chart pocket
- Carbon steel frame provides strength and durability
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Hammertone paint provides a long lasting and durable finish
- Comfortable padded armrests
- Maximum Weight Capacity 250 lb – EVENLY DISTRIBUTED
- Meets Medicare code K0001

Warranty

One (1) year limited warranty



**This product is not made with natural rubber latex*



Fixed Full Arm,
Swingaway Footrest



Detachable Desk Arm,
Swingaway Footrest



Detachable Desk Arm,
Elevating Legrest

Advantage		Dimensions and Weights									
		Seat			Back	Overall				Wt.	Max Wt.
											
3H010100	Fixed Full Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010110	Fixed Full Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010120	Detachable Desk Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010130	Detachable Desk Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010140	Detachable Full Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010150	Detachable Full Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010200	Fixed Full Arm, Swingaway Footrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010220	Detachable Desk Arm, Swingaway Footrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010230	Detachable Desk Arm, Elevating Legrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010240	Detachable Full Arm, Swingaway Footrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010250	Detachable Full Arm, Elevating Legrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H010320	Detachable Desk Arm, Swingaway Footrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	250 lb
3H010330	Detachable Desk Arm, Elevating Legrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	250 lb
3H010340	Detachable Full Arm, Swingaway Footrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	250 lb
3H010350	Detachable Full Arm, Elevating Legrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	250 lb
3H011100	Vinyl, Fixed Full Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H011110	Vinyl, Fixed Full Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H011120	Vinyl, Detachable Desk Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H011130	Vinyl, Detachable Desk Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H011220	Vinyl, Detachable Desk Arm, Swingaway Footrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H011230	Vinyl, Detachable Desk Arm, Elevating Legrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	250 lb
3H011320	Vinyl, Detachable Desk Arm, Swingaway Footrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	250 lb
3H011330	Vinyl, Detachable Desk Arm, Elevating Legrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	250 lb

Advantage Series – Common Accessories & Parts		
90763064	Heel loop, pair	E0951
90763062	Wheel lock extension, pair	E0961
B10320F	O2-IV holder	E2208
B10325	O2 cylinder holder	E2208
GF1844W	Chrome-plated IV pole for W/C 1 hook	K0105
EJO2BAG07	O2 holder bag	E2208
90763065	Anti-tippers, pair 20" seat height	E0971
90763066	Anti-tippers, pair 18" seat height	E0971
92523850	Anti-fold bar (chromed) lower frame	
9B761770	Anti-theft cap for legrest for Traveler chairs with 1.5" spacing	