

B700 BARIATRIC BED



Full-Electric B700 Bariatric Bed Shown with Head Half Rails.

B700 Bariatric Bed

- Height travel range: 17" to 27"
- Safe work load: 700 pounds*
- Full-electric, three-function operation
- DC Electronics that are quieter and offer a faster Hi-Lo function than traditional AC
- DC Electronics are splashproof and offer a lower voltage making them more reliable and less prone to electronic failures
- Designed to meet the requirements of applicable HCPCS codes
- Split-frame with removable ends for easy transport
- High-strength, steel-reinforced, truss frame
- Heavy-duty grid sleep deck
- Durable powder-coated finish
- Simple assembly: no special tools required
- Self-wrenching nut: assembled for easy part replacement
- Ball-screw bed ends

* 700 lbs. evenly distributed total weight capacity including the patient, mattress, rails and accessories

- Overall bed dimensions: 46.5" width (with rails) and 87" length
- Sleeping surface dimensions: 42" width and 80" length
- Total Bed Weight: 295 lbs
- Optional Swing-down style trigger release rail requires no tool assembly and designed with quick release for easy removal

Warranty:

- Two year warranty for defects in workmanship and materials of mechanical components, frame, and electronics
- During the warranty period, defective items will be repaired or replaced at manufacturer's option at no charge



5" locking caster

ITEM	Descriptions	HCPCS Code
ABL-B700	Full-Electric Bed Only	E0301
ABL-B700-PKG	Bed, Head Half Rails, Bariatric Mattress	E0303
650-1501-1633	Bariatric Mattress	E0184
A8340	Head Half Rails, pr	E0305
A8510	Foot Half Rails, pr	E0305



1.678.291.3200 Fax 1.770.368.4954

Information contained herein is subject to change.

The most current and complete product information can be found on our website

www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved.

Basic American Medical Products is a trademark of GF Health Products, Inc.

GF Health Products, Inc. is an ISO 13485 Certified Company

GF1200043RevA12