

Everest & Jennings®

WHEELCHAIRS

Setting standards in quality
Enhanced transfer and comfort
Manufactured for durability





Everest & Jennings®

Value:

Graham-Field® is known for a broad product offering at affordable pricing; without sacrifice to quality. Like you, we understand that the cost of a product is more than a price tag. The cost also includes the time and effort required to sell and service your customers. At Graham-Field we are proud to offer the value found in the Everest & Jennings™ line of wheelchairs.

Quality:

E&J™ has maintained a commitment to provide our customers quality manual folding wheelchairs since 1933. In a time when many manufacturers are producing products at the lowest possible cost, often cutting corners on quality, E&J is standing by its dedication to quality through product design, specification control, testing and inspection. It is easy to see why the E&J Traveler® name is synonymous with quality.

Service:

Our products are often only a day away with centralized distribution located in the metro areas of Atlanta, Georgia; St. Louis, Missouri; and Los Angeles, California. E&J also carries a full line of parts and upholstery to meet your service needs. In addition, our customer service department and technical service teams are always ready to help with your product needs.



Manual Folding Wheelchairs

Advantage®

The Advantage wheelchair provides comfort and durability at an affordable price. Designed as a low maintenance wheelchair.

Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth
- Black padded nylon and vinyl upholstery available with chart pocket
- Carbon steel frame provides strength and durability
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Hammertone paint provides a long lasting and durable finish
- Comfortable padded armrests
- Maximum Weight Capacity – 250 lb – EVENLY DISTRIBUTED
- Meets Medicare code K0001
- One (1) year limited warranty



Traveler® SE

The Traveler SE features a hammertone frame with saddle-joint tubular construction with dual axle for an adjustable 2" seat height and tool-free, industry compatible, front rigging.

Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth and 20" by 18" seat depth
- Black vinyl upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Two-position axle and caster fork provide 2" seat height adjustment
- Maximum Weight Capacity 250 lb & 300 lb on 20"x18" models - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002
- One (1) year limited warranty



Traveler® LX

The Traveler LX incorporates the legendary Everest & Jennings wheelchair chrome finish frame with saddle-joint tubular construction with dual axle for an adjustable 2" seat height and tool-free industry compatible front rigging. The premium vinyl upholstery is reinforced with an extra foam and heavy duty inner liner for enhanced comfort and durability.

Standard Features

- Available in 16", 18" and 20" seat widths
- Blue premium vinyl upholstery with standard chart pocket on back
- Chrome finish provides a long lasting and easy to clean finish
- Tool-free front rigging with heel loops standard on swingaway footrests
- Gray non-marking tires
- Two-position axle and caster fork provide 2" seat height and adjustment
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002
- Two (2) year limited warranty



Bariatric Folding Wheelchairs

Traveler® HD

The Traveler HD, featuring a heavy-duty reinforced frame, provides a solution for users weighing up to 500 lb.

The Traveler HD comes with double crossbraces, heavy-duty 5/8" diameter axles, and a long lasting, embossed back and seat upholstery with a double inner liner for durability.

Standard Features

- Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard)
- Embossed black vinyl upholstery with double inner lining for durability
- Hammertone paint provides a long lasting and durable finish
- Reinforced frame at all critical load-bearing points
- Chart pocket standard on upholstery
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



Paramount XD

The Paramount XD is ideal for the bariatric wheelchair user. It features 26" and 30" seat widths as well as a weight capacity of 650 lb. The Paramount XD is constructed of 1" tubing to meet high weight capacity demands.

Standard Features

- Available in 26" and 30" seat widths by 19 1/2" seat depth
- Swingaway footrest and elevating legrest options
- Maximum weight capacity of 650 lb - EVENLY DISTRIBUTED
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Reinforced welds and gussets on frame
- 8" x 2" Solid casters
- 24" Solid rear wheels
- Meets Medicare code K0007
- Five (5) year limited warranty



Traveler® HTC

The Traveler HTC (Hospital Transport Chair) provides an alternative solution to the hospital/clinical transport need by incorporating special attendant and inventory control features with a heavy duty frame reinforced at all critical load-bearing points.

The Traveler HTC comes with an attendant push handle and rear mounted brakes that offer easier transport for the caregiver. In addition, the full width push handle acts as an anti-fold bar.

Standard Features

- Available in 20", 22" and 24" seat widths, 18" seat depth (1" longer than standard)
- Fixed attendant push handle also functions as an anti-fold bar
- Includes a fixed anti-theft post ideal in hospital/clinical settings
- Rear mounted attendant wheel locks
- Fixed desk arms and front rigging lock nuts ensure inventory control
- Anti-tippers standard
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Black vinyl upholstery with chart pocket on back
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007
- Five (5) year limited warranty



Specialty Folding Wheelchairs

Advantage® Recliner

The Advantage Recliner comes with a higher back and removable full support headrest. The headrest includes an additional neck support that adjusts easily to fit every need. The durable hammertone frame is designed to support up to 300 lb users featuring a full reclining back range from 90-170 degrees.

Standard Features

- Available in 16", 18" and 20" seat widths, with 17" seat depth
- Manual full reclining back
- Removable full support headrest
- Detachable desk or full arms
- Hammertone finish provides durability and style
- Embossed black vinyl upholstery
- Anti-tippers standard
- Comes with spreader back bar to support upholstery and works as anti-fold device
- Industry compatible tool-free adjustable elevating front rigging
- Chart pocket standard on back
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001/ E1226 (Recliner Back) / E0966 (Head Rest Extension)
- Five (5) year limited warranty



Rehab Shower Commode

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or used as a portable, self-contained commode.

Standard Features

- Available in 14 1/2" and 17" High back and Low back, caster pin locks and 5" rear locking caster models
- Completely removable full length arms
- Swingaway, detachable footrests
- Non-corrosive composite footplates
- Recessed seat and frame side contours
- Removable, hygienic plastic pail and lid standard
- Oval opening padded seat comes as standard
- Cushioned vinyl seat is sealed against moisture penetration
- Maximum weight capacity 250 lb – EVENLY DISTRIBUTED
- Meets Medicare code E0240
- Two (2) year limited warranty



1.770.368.4700

Information contained herein is subject to change.
The most current and complete product information can be found on our website.
www.grahamfield.com

Lightweight Folding Wheelchairs

Traveler® L3

The Traveler L3 lightweight frame weighs less than 36 lb making it the ideal wheelchair built for short and long term use. Designed with dual axles and multi-position caster forks, this wheelchair provides an adjustable seat height of 18" or 20". The wheelchair comes standard with non-marking polyurethane tires and industry compatible front rigging.



Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth
- Black nylon upholstery with standard chart pocket on back
- Non-marking polyurethane rear tires
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Two-position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003
- Five (5) year limited warranty

Traveler® L4

The Traveler L4 is a high-strength lightweight K0004 Medicare coded wheelchair. The Traveler L4 has dual axle and multi-position caster forks for an adjustable seat height of 18 or 20 inches. The wheelchair comes standard with non-marking polyurethane tires, anti-tippers and a chart pocket.



Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth. Also 18" and 20" seat widths with 18 seat depth
- Quick release axle is available for some models
- Black nylon upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Anti-tippers come standard
- Flip-back desk and full length arms to ease lateral transfer
- Height adjustable desk arms on some models
- Industry compatible tool-free front rigging with standard heel loops on swingaway footrests
- Two-position axle and caster fork provides 2" seat height adjustment
- Maximum weight capacity 250 lb EVENLY DISTRIBUTED (300 lb on models with 18" seat depth)
- Meets Medicare code K0004
- Five (5) year limited warranty

E&J Navigator

The E&J Navigator easily converts from a classic transport chair to a lightweight self-propelled wheelchair, using the included 24" quick-release axle rear wheels. The rear wheels are easily attached or removed with the push of a button. No additional tools are required.

Other features include convenient flip back desk arms and swingaway footrests for easy patient transfer. The back of the Navigator folds down making transport and storage easy. The compact size, flexible configurations, lightweight and high strength frame make this a transport and lightweight wheelchair all in one.

Standard Features

- 18" seat width
- Transport chair and lightweight wheelchair all in one
- Available in 3 stylish colors – Black, Blue and Red
- Aluminum frame
- Black padded nylon upholstery
- Convenient flip-back desk arms
- Fold down back makes transport and storage easy
- Non marking polyurethane tires
- Weighs 33.5 lb with 24" rear wheels installed, 23 lb without 24" rear wheels (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003
- Standard rear mount anti-tippers
- Two (2) year limited warranty



Transport Chairs

Aluminum Transport

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device and known for its durability. Its compact size and lighter weight allow a companion to lift and store the chair with ease.



Rear mounted attendant locks make it easy for a companion to engage.

Standard Features

- 17" and 19" seat widths
- Lightweight aluminum frame available in stylish colors options of Aqua, Black, Red, Blue, Silver and Pink
- Fold down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Footrests are industry compatible and offer easy tool-free adjustment
- Full length padded arms provide user comfort
- Weighs 20 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code E10385
- One (1) year limited warranty

Steel Transport

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device and known for its durability. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

Rear mounted attendant locks make it easy for a companion to engage.

Standard Features

- 17" and 19" seat widths
- Steel frame is finished in a durable hammertone coating
- Fold down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Footrests are industry compatible and offer easy tool-free adjustment
- Full length padded arms provide user comfort
- Weighs 24 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code E10385
- One (1) year limited warranty



Bariatric Transport

The Everest & Jennings Bariatric Transport Chair is designed as an affordable bariatric patient transfer device. With its 24" width and 400 lb weight capacity, this chair satisfies additional compact transport needs with extra capacity. Handle mounted attendant locks make it easy for a companion to engage.

Standard Features

- 24" seat width
- Lightweight aluminum frame with silver epoxy painted silver color finish
- Fold down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Standard Carry straps on seat allow for ease of lifting and folding
- Padded flip-back desk arms provide user comfort
- Weighs 31 lb (w/o front rigging)
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Meets Medicare code E1039
- Three (3) year limited warranty



Accessories

Anti-tippers	
	
90763565	Anti-tipper, Advantage
90763065	Anti-tipper, Traveler SE & HD
90763065	Anti-tipper, Traveler L3

IV/O2 Holders	
	
B10320F	O2/IV Holder Combination (Advantage, SE, L3, L4, HD, LX, Paramount XD)
B10325F	O2 Holder (Advantage, SE, L3, L4, HD, LX, Paramount XD)
GF1844W	IV Holder (Advantage, SE, L3 and Traveler Chromed Recliner)

Positioning Belts	
	
92133040	Positioning Belt, Metal Buckle
300-505	Positioning Belt, Auto Style
EJ76-BELT	Positioning Belt, Plastic Buckle

O2 Cylinder Holder Bag	
	
EJ02BAG07	O2 Bag

Footrest	
	
90763030	Footrest with Composite Footplate
90763430	Footrest with Aluminum Footplate

High Pivot Legrest	
	
90763435	High Pivot Elevating Legrest with Aluminum Footplate, Hammertone Finish
90763135	High Pivot Elevating Legrest with Aluminum Footplate, Chrome Finish

Legrest	
	
90763031	Elevating Legrest with Composite Footplate
90763431	Elevating Legrest with Aluminum Footplate

Padded Seat	
	
90635070	Oval Cutout
90635175	Open Front (C-shaped)

Item No.	Description
Advantage (K0001)	
3H010220 / 3H010230	Advantage 16X16 Desk Arms, Footrest / Legrest
3H010120 / 3H010130	Advantage 18X16 Desk Arms, Footrest / Legrest
3H010320 / 3H010330	Advantage 20X16 Desk Arms, Footrest / Legrest
Traveler SE (K0001/K0002)	
3E010220 / 3E010230	Traveler SE 16X16 Desk Arms, Footrest / Legrest
3E010120 / 3E010130	Traveler SE 18X16 Desk Arms, Footrest / Legrest
3E010320 / 3E010330	Traveler SE 20X16 Desk Arms, Footrest / Legrest
Traveler LX (K0001/K0002)	
3E020220 / 3E020230	Traveler LX 16X16 Desk Arms, Footrest / Legrest
3E020120 / 3E020130	Traveler LX 18X16 Desk Arms, Footrest / Legrest
3E020320 / 3E020330	Traveler LX 20X16 Desk Arms, Footrest / Legrest
Traveler L3 (K0003)	
3F010220 / 3F010230	Traveler L3 16X16 Desk Arms, Footrest / Legrest
3F010120 / 3F010130	Traveler L3 18X16 Desk Arms, Footrest / Legrest
3F010320 / 3F010330	Traveler L3 20X16 Desk Arms, Footrest / Legrest
Traveler L4 (K0004)	
3F020220 / 3F020230	Traveler L4 16X16 Desk Arms, Footrest / Legrest
3F020120 / 3F020130	Traveler L4 18X16 Desk Arms, Footrest / Legrest
3F020320 / 3F020330	Traveler L4 20X16 Desk Arms, Footrest / Legrest
Navigator Wheelchair & Transport Chair	
EJNAV-18BK / BL/ RD	E&J NAVIGATOR 18W BLACK / Blue / Red
Traveler HD Heavy Duty (K0007)	
3G010320 / 3G010330	Traveler HD 20X18 Desk Arms, Footrest / Legrest
3G010420 / 3G010430	Traveler HD 22X18 Desk Arms, Footrest / Legrest
3G010520 / 3G010530	Traveler HD 24X18 Desk Arms, Footrest / Legrest
Traveler Hospital Transport Chair (HTC)	
3J010330	Traveler HTC 20X18 Desk Arms, Elevating Legrest
3J010430	Traveler HTC 22X18 Desk Arms, Elevating Legrest
3J010530	Traveler HTC 24X18 Desk Arms, Elevating Legrest
Paramount XD Heavy Duty (K0007)	
5PX10620 / 5PX10630	Paramount XD 26W Desk Arms, Footrest / Legrest
5PX10820 / 5PX10830	Paramount XD 30W Desk Arms, Footrest / Legrest
Aluminum Transport Chairs	
EJ789-1 / EJ780-1	Transport Chair Green 17" / 19"
EJ781-1 / EJ782-1	Transport Chair Red 17" / 19"
EJ783-1 / EJ784-1	Transport Chair Silver 17" / 19"
EJ785-1 / EJ786-1	Transport Chair Black 17" / 19"
EJ787-1 / EJ788-1	Transport Chair Blue 17" / 19"
EJ790-1 / EJ791-1	Transport Chair Pink 17" / 19"
EJ793-1	Transport Chair Aqua 19"
Steel Transport Chair	
EJ796-1	Transport Chair Silver Vein Steel 19"
Bariatric Transport	
EJ777-3	Transport Bariatric Chair 24W
Advantage Reclining Wheelchair	
3K010230 / 3K010250	Advantage Recliner 16X17 Legrest, Desk Arms / Full Arms
3K010130 / 3K010150	Advantage Recliner 18X17 Legrest, Desk Arms / Full Arms
3K010330 / 3K010350	Advantage Recliner 20X17 Legrest, Desk Arms / Full Arms
Rehab Shower Commode	
12022010 / 12022011	14 1/2" LoBack / 14 1/2" LoBack, Caster Pin Locks
12022013	14 1/2" LoBack with 5" Rear Locking Caster
12022110 / 12022111	17 1/2" HiBack / 17 1/2" HiBack, Caster Pin Locks

*For a complete list of sizes and options available, please visit www.grahamfield.com