

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or used as a portable, self-contained commode.



Optional C-shaped seat available



Standard Features

- Completely removable full length arms for ease of transfer
- Swingaway, detachable footrests permit close approaches to beds, toilets and wheelchairs
- Recessed seat and frame side contours allow easy access for toileting activities while seated
- Available with durable 24" one piece composite wheels with integral handrims and solid rubber tires or 5" rear casters
- Removable, hygienic plastic pail and lid are included
- Cushioned vinyl seat is sealed against moisture penetration and vented for optimal foam compression
- Non-corrosive composite footplates
- Maximum weight capacity 250 lb **EVENLY DISTRIBUTED** on models with 14 1/4" & 17 1/4" back height
- Maximum weight capacity 300 lb **EVENLY DISTRIBUTED** on models with 16" & 19" back height
- Chair comes standard with oval opening seat
C-shaped seat is an optional item
- Meets Medicare code E0240



Warranty

Two (2) year limited warranty





Rehab Shower Commode		Dimensions and Weights								
		Seat			Back	Overall			Weight	
		Width	Depth	Height	Height	Depth	Height	Width	Without front rigging	Max Capacity
12022010	14 1/4" LoBack	17"	18"	23"	14 1/4"	41 1/2"	37 1/2"	24 1/2"	43 lb	250 lb
12022011	14 1/4" LoBack, Caster Pin Locks	17"	18"	23"	14 1/4"	41 1/2"	37 1/2"	24 1/2"	44 lb	250 lb
12022013	14 1/4" LoBack with 5" Rear Locking Caster	17"	18"	23"	14 1/4"	41 1/2"	37 1/2"	24 1/2"	36 lb	250 lb
12022110	17 1/4" HiBack	17"	18"	23"	17 1/4"	41 1/2"	40 1/2"	24 1/2"	43 lb	250 lb
12022111	17 1/4" HiBack, Caster Pin Locks	17"	18"	23"	17 1/4"	41 1/2"	40 1/2"	24 1/2"	44 lb	250 lb
12023010	16" LoBack	20"	19"	23 1/2"	16"	42 3/4"	40"	28 1/4"	48 lb	300 lb
12023110	19" HiBack	20"	19"	23 1/2"	19"	42 3/4"	43"	28 1/4"	50 lb	300 lb

All dimensions are ±1/4"
*Without Front Rigging

Rehab Shower Commode – Common Accessories & Parts		HCPSC Code
90835070	Oval cutout seat for 17 1/4" and 14 1/4"	
90835071	Oval cutout seat for 16"	
90835175	Open front (C-shaped) seat for 17 1/4" and 14 1/4"	
90H43170	Commode pail & lid	
90761752	Footrest kit stainless steel (pair)	K0045
90761752L	Footrest kit stainless steel (pair) (long)	K0045