




Advantage[®]

PARTS CATALOG

3H01 MODELS



Everest & Jennings[®]
by graham-field

CONTENTS

SIDEFRA^{ME} (FIXED ARM), CROSSBRACE, ARM PAD KIT3

SIDEFRA^{ME} (DETACHABLE ARM), CROSSBRACE.....4

ARM, DETACHABLE5

CASTER.....5

REAR WHEEL5

FRONT RIGGING, DETACHABLE6

UPHOLSTERY7


ACCESSORIES7

TO ORDER

Tel: 1.770.368.4700 • Fax: 1.770.368.2386

Returns: Parts may not be returned without prior authorization from GF Health Products, Inc. Custom parts cannot be returned for credit. Please refer to GF Health Products, Inc. [Terms and Conditions](#) for further information.

SIDEFRADE (FIXED ARM), CROSSBRACE, ARM PAD KIT



1 Sideframe, fixed arm: Includes side panel

Right 90763501
Left 90763502

2 Crossbrace, pair: Includes center bolt and nut


16" W 907635H1
18" W 907635H2

3 Bottom inner rail kit: Includes bottom inner rail (2 each) and screw and nut that hold bottom inner rail (4 each)..... 90763503

4 Grip (10 each) 90763005


5 Seat rail guide (4 each) 90763004

6 Tipping sleeve (10 each)..... 90763006



Arm pad kit: Includes 2 black
full-length arm pads and
mounting screw..... 90763524

SIDEFRADE (DETACHABLE ARM), CROSSBRACE



1 Sideframe, detachable arm

Right 90763511
Left 90763512

2 Crossbrace, pair: Includes center bolt and nut

16" W 907635C1
18" W 907635C2
20" W 907635C3

3 Bottom inner rail kit: Includes bottom inner rail (2 each); screw and nut that hold bottom inner rail (4 each); and sleeve (2 each)

..... 90763503

4 Seat slide kit:

Single telescopic: Includes seat slide (2 each) and bushing (2 each)

16" and 18" W 90763513

Double telescopic: Includes seat slide (2 each)

20" W 90763517


5 Grip (10 each) 90763005

6 Seat rail guide, rear (2 each) 90763015

7 Tipping sleeve (10 each) 90763006

8 Inner seat rail (2 each) 90763514

ARM, DETACHABLE



- 1 Arm, detachable fixed-height, pair:**
Includes right and left arm with side panel and spring lock (does not include arm pad)


Desk length..... 90763020
Full length..... 90763021

- 2 Arm pad, pair:** Includes two black arm pads and mounting screws

Desk length..... 90763022
Full length..... 90763023

- 3 Spring lock (10 each)** 90763024


CASTER



Complete caster (fork and wheel),
8" x 1" solid rubber tire..... 90763540

- 1 Caster wheel, axle and nut** 90763541
2 Caster fork, stem nut and dust cap 90763542
3 Caster stem bearing (10 each) 90763543
4 Caster dust cap (10 each)..... 90763544
5 Caster axle and nut (10 each)..... 90763545

REAR WHEEL



- 1 Rear wheel with handrim, 24" x 1",**
solid rubber tire: Includes hub cap and bearing; does not include axle or nut 90763550


- 2 Hub cap (10 each)**..... 90763051

- 3 Wheel bearing (10 each)** 90763552

- 4 Rear axle and nut (10 each)**


Fixed-arm sideframe 90763553
Detachable-arm sideframe 90763554

FRONT RIGGING, DETACHABLE




Footrest, pair, all widths
Includes right and left hanger with release
and composite footplate with extension
and heel loop 90763030

Composite footplate, pair
(Fits all front rigging) 90763033




Low-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
composite footplate with extension, and
legrest panel (does not include
heel loop) 90763031

Elevating legrest panel kit
(Fits all elevating legrests)
Includes 2 black padded panels
and hardware 90763032



High-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
aluminum footplate with extension, and
legrest panel (does not include
heel loop) 90763435




NYLON UPHOLSTERY

Back upholstery, black, with built-in chart holder: Includes hardware

- Nylon, fixed arm, 16" wide 907635D1
- Nylon, fixed arm, 18" wide 907635D2
- Nylon, detachable arm, 16" wide ... 907635A1
- Nylon, detachable arm, 18" wide ... 907635A2
- Nylon, detachable arm, 20" wide ... 907635A3

Seat upholstery, black: Includes mounting strips and hardware

- Nylon, 16" wide 907635B1
- Nylon, 18" wide 907635B2
- Nylon, 20" wide 907635B3



VINYL UPHOLSTERY

Back upholstery, black, with built-in chart holder: Includes hardware

- Vinyl, fixed arm, 16" wide 907635F1
- Vinyl, fixed arm, 18" wide 907635F2
- Vinyl, detachable arm, 16" wide 907635E1
- Vinyl, detachable arm, 18" wide 907635E2
- Vinyl, detachable arm, 20" wide 907635E3

Seat upholstery, black: Includes mounting strips and hardware

- Vinyl, 16" wide 907635G1
- Vinyl, 18" wide 907635G2
- Vinyl, 20" wide 907635G3

ACCESSORIES

Advantage accessories (not shown)

- 1 Heel loop, pair 90763064
- 2 Wheel lock, pair: Includes rubber tips
 - Fixed-arm sideframe**
 - Push to lock 90763070
 - Pull to lock 90763071
 - Detachable-arm sideframe**
 - Push to lock 90763060
 - Pull to lock 90763061
- 3 Rubber tip, wheel lock (10 each) 90763063

- 4 Wheel lock extension, pair 90763062
- 5 Positioning belt (metal buckle) 92133040
- 6 Anti-Fold bar 92523850
- 7 O2-IV holder B10320F
- 8 O2 cylinder holder B10325
- 9 Anti-Tippers, pair:
 - Fixed-arm sideframe 90763566
 - Detachable-arm sideframe 90763565



1.770.368.4700

Information contained herein is subject to change.
The most current and complete product information can be found on our website.
www.grahamfield.com

Follow us on



© 2009, GF Health Products, Inc. All Rights Reserved.
Advantage and Everest & Jennings are trademarks of GF Health Products, Inc.,
One Graham-Field Way, Atlanta GA 30340-3140

GF Health Products, Inc. is an ISO 13485 Certified Company.