

OVERBED TABLES



The Lumex Overbed Table combines an attractive, ergonomic design with features that reduce infection control issues by minimizing areas in which dust and dirt can accumulate. Manufactured with the lowest base in the industry, our Overbed tables move easily under any hospital bed or recliner.



Features

- Lowest profile base on the market designed to fit under any hospital bed
- Exposed steel base with no plastic cover allows for easy cleaning of base and casters
- Laminate or Thermofoil top that provides durability
- Available in 4 attractive finishes: Medium Oak, Medium Cherry, Hardrock Maple and Candlelight
- Seamless Spill guard top holds up to 36 oz
- Easy, one-hand activated secondary surface reduces infection control issues and eliminates hidden spaces between surfaces for bacteria to collect
- Integrated cup holders can hold soda cans, water bottles and medication
- Optional Vanity drawer with mirror slides out from either side of the table
- Integrated accessory rail on both sides of column to easily transfer accessories such as cord wraps, IV hooks, urinal hooks and storage baskets from headwall rails
- Shock-absorbent gas spring cylinder for easy height adjustment
- Designed to support up to 100 lb of evenly distributed weight
- Two (2) Year Limited Warranty

| GF Stock Code | Description |
|---------------|--|
| GF9000-LSTMC | OVERBED TABLE - SINGLE TOP, LAMINATE, MEDIUM CHERRY GRAPHITE GREY U-BASE |
| GF9000-LSTHM | OVERBED TABLE - SINGLE TOP, LAMINATE, HARDROCK MAPLE, GRAPHITE GREY U-BASE |
| GF9000-LSTCL | OVERBED TABLE - SINGLE TOP, LAMINATE, CANDELLIGHT, GRAPHITE GREY U-BASE |
| GF9000-TSTM0 | OVERBED TABLE - SINGLE TOP, THERMOFOIL, MEDIUM OAK, GRAPHITE GREY U-BASE |
| GF9000-TSTM C | OVERBED TABLE - SINGLE TOP, THERMOFOIL, MEDIUM CHERRY GRAPHITE GREY U-BASE |
| GF9000-TSTHM | OVERBED TABLE - SINGLE TOP, THERMOFOIL, HARDROCK MAPLE, GRAPHITE GREY U-BASE |
| GF9000-TSTCL | OVERBED TABLE - SINGLE TOP, THERMOFOIL, CANDELLIGHT, GRAPHITE GREY U-BASE |
| GF9000-TDTM0 | OVERBED TABLE - FLIP TOP, THERMOFOIL, MEDIUM OAK, GRAPHITE GREY U-BASE |
| GF9000-TDTM C | OVERBED TABLE - FLIP TOP, THERMOFOIL, MEDIUM CHERRY GRAPHITE GREY U-BASE |
| GF9000-TDTHM | OVERBED TABLE - FLIP TOP, THERMOFOIL, HARDROCK MAPLE, GRAPHITE GREY U-BASE |
| GF9000-TDTCL | OVERBED TABLE - FLIP TOP, THERMOFOIL, CANDELLIGHT, GRAPHITE GREY U-BASE |
| GF9000-BSK | OVERBED TABLE - STORAGE BASKET - 10x6x9 WITH ADAPTER |
| GF9000-VNTY | OVERBED TABLE - VANITY DRAWER WITH MIRROR |
| GF9000-UH | OVERBED TABLE - URINAL HOOK WITH ADAPTER |
| GF9000-CST | OVERBED TABLE - NON-LOCKING CASTER |

Standard Top Finishes



Accessories



Standard Dimensions:

Single Table Top: 32" L x 17" W (81.3 cm x 43.2 cm)
 Flip Table Top: 13.5" L x 17" W (34.3 cm x 43.2 cm)
 Storage Table Top: 41.5" L x 19" W (105.4 cm x 48.3 cm)
 Height Range: 28" Low - 44" High (71 cm - 111.8 cm)

1.770.368.4700

Information contained herein is subject to change.

Contact your Graham-Field Representative at

1-800-347-5678 or visit our website.

www.grahamfield.com

Follow us on



© 2019, GF Health Products, Inc. All Rights Reserved.
 GF Health Products, Inc. is an ISO 13485 Certified Company.
 GF1900018RevA19