

Multipurpose tray is ideal for eating, reading, and writing while seated in wheelchair or resting in bed



Standard Features:

- Ergonomic design with contoured cutout and comfortable beanbag style base
- Easy to clean composite top measures 15.75" x 13" x 2" including 2.5" diameter (40 cm x 33 cm x 5 cm including 6.5 cm diameter) engraved circumferences for sitting cups or bottles
- Washable nylon cover includes interior 6" (15.24 cm) zippered small pocket on top and side pocket
- Convenient carry handle attached
- Composite tray can be cleaned with soap and water or diluted bleach
- Cover can be removed and hand washed with soap and water, re-attached to fastener tabs
- Total weight: 1.3 lb (600g)
- Weight capacity: 10 lb (4.53kg) evenly distributed



Warranty:

1 (one) year limited warranty

Item	UOM
EJTRAY18	1 ea



**This product is not made with natural rubber latex*