

The Lumex Knee Walker is ideal for individuals recovering from sprains, foot surgery, amputation or any foot or ankle injury where the user cannot bear weight. It can also be a more comfortable alternative to crutches which with extended use can lead to underarm and/or back pain. The Lumex Knee Walker can help with a speedy recovery and help an individual maintain a more independent and active lifestyle.



Features:

- Ideal for individuals recovering from foot surgery, Achilles injuries, breaks and/or sprains of the foot
- Knee walker turning system can be steered for increased maneuverability
- Built-in shock absorber is designed to absorb shock impulses during movement making it a smoother ride for the user
- Deluxe, padded, contoured polyurethane knee pad is designed for ultimate comfort
- Braking system allows for quick and extended stopping
- Easy, tool-free adjustable handle and knee pad accommodate various user heights
- Robust 8" non-marking polyurethane casters are ideal for indoor/outdoor use
- Durable, aluminum frame with attractive color options of Apple Red, Black, and Purple
- Convenient, removable front storage pouch is ideal for storage of personal items
- 300 lb (136 kg) maximum weight capacity – EVENLY DISTRIBUTED
- Lifetime (3 Years) Limited Warranty
- Designed to meet the requirements of E0118



Item	Description	UOM
LX8000A	Apple Red	1 ea
LX8000BK	Black	1 ea
LX8000PU	Purple	1 ea

Specifications

Overall Width.....	19"
Overall Length.....	31"
Minimum Handle Height.....	34.3"
Maximum Handle Height.....	41.3"
Knee Pad Dimensions.....	14" x 7" x 2.4" (LxWxH)
Knee Pad Minimum Height.....	17.75"
Knee Pad Maximum Height.....	22.75"
Maximum Weight Capacity.....	300 lb (136 kg) EVENLY DISTRIBUTED
Product Weight.....	17.75 lb (8.05 kg)



Braking System



Shock Absorber