



## CLINICAL CARE RECLINERS

# CLINICAL CARE RECLINER - WIDE WITH DROP ARM

Three positions for patient comfort. The fourth position is a caregiver-activated Trendelenburg (or Shock) Position.

Two easy to use Self Storing Drop Arms allow for a safer lateral transfer for both the patient and caregiver

Two self storing Side Tables for patient prep/convenience; smooth surface allows for easy cleaning

Dual gas cylinders give smooth power assistance during Trendelenburg activation.

Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster.

Wall saver Armrest Tray increases durability of the armrests and protect walls from marring (Recliner is also available with Fully Upholstered Arms)

Side Panels are easily removed for cleaning and maintenance

All Recliner components can be replaced greatly enhancing the lifespan of the recliner and reducing its cost in use over time

Wide selection of fabrics and colors. See the [Swatch Card](#) for details.



### Options:

- Heat and Massage
- Heat Only feature
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow



Multiple Positions for Enhanced Patient Comfort

Model (xxxx represents fabric colors)	Wall Saver Armrest Edges	Fully Upholstered Arms	Heat & Massage	Heat Only	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Width between side tables in upright position		Depth fully reclined (Trendelenburg)		Weight Capacity EVENLY DISTRIBUTED		UOM
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
FR587WD	●				33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WDH	●		●		33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WDHO	●				33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WDxxxxPA	●				33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WDHxxxxPA	●		●		33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WDHOxxxxPA	●		●		33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2018, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.