


DIAGNOSTIC INSTRUMENTS

ORTHOPEDIC GONIOMETER

Features

- The orthopedic goniometer is used to measure a joint's range of motion (ROM). There are two "arms"—one that is stationary and one that is movable—that are hinged together
- Plastic



| | Length | | UOM |
|---------|--------|------|------|
| | inch | cm | |
| 12-1000 | 12 | 30.5 | 1 ea |
| 12-1001 | 8 | 20.3 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.

graham
field

GF2100150RevC22