

## EXAM ROOM SUPPLIES

# STANDARD MAGNIFYING BINOCULAR LOUPE

### Features

- Single lens, with no center post, increases viewing area for clear, comfortable, unobstructed vision
- Shatter-resistant
- Magnification brings smallest details up large and sharp
- Interchangeable lens plate units snap on and off instantly without screws or tools
- Convenient flip-up, flip-down feature for regular vision and magnification
- Can be worn over regular or safety
- Lens is hard-coated acrylic
- Lightweight: 4 ½ oz (127.6 g)
- Self locking, vinyl-lined, adjustable headband



	Power	Working Distance		UOM
		inch	cm	
1200	1.75X	14	35.6	1 ea
1201	2.25X	8	20.3	1 ea
1202	2.75X	6	15.2	1 ea
1203	3.5X	4	10.2	1 ea



GF Health Products, Inc. is an  
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2005, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF2100199Rev/C22