

LABORATORY

# LENS PAPER

Lens paper is a soft, lint-free tissue paper, and used to clean the lenses of microscopes.

**Features**

- Soft, thin, silky paper sheets
- (50) 4" x 6" (10.2 x 15.2 cm) sheets per book

**Technical specifications:**

- AA-grade wood pulp
- 0.02mm thick
- 13gm/square meter



	overall Width		overall Length		UOM
	inch	cm	inch	cm	
3719	4	10.2	6	15.2	12 books/pk



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2019, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.