

## MOBILITY ACCESSORIES

# PLATFORM WALKER ATTACHMENT

### Features

- Designed for users with limited strength or dexterity who are unable to grip the walker
- Adjustable strap comfortably secures arm in the forearm rest
- Padded forearm rest can be pivoted and angled to fit user comfortably and to assist in lifting walker
- Padded handgrip provides comfort
- Universal clamp fits most walkers
- Tool-free installation and adjustment allows quick fitting for most patient sizes and types
- Designed to meet the requirements of HCPCS code: E0154
- Three (3) year limited warranty



*Platform Attachment is recommended for use on Lumex Walkers. (We cannot guarantee the performance of this product on other industry walkers.)*

		UOM
6132B	Platform Walker Attachment	1 pr



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2018, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

\* It is recommended that the platform walker attachments be installed on both sides for better stability.

**graham**  
**field**

GF-1800-063RevC22