

## OVERBED TABLES

Manufactured with the lowest base in the industry, our Overbed tables move easily under most hospital beds or recliners.

### Features

- Lowest profile base on the market
- Exposed steel base with no plastic cover allows for easy cleaning of base and casters
- Laminate or Thermofoil top that provides durability
- Available in 4 attractive finishes: Candlelight, Hardrock Maple, Medium Cherry, and Medium Oak
- Easy, one-hand activated secondary surface reduces infection control issues and eliminates hidden spaces between surfaces for bacteria to collect
- Integrated cup holders
- Optional Vanity drawer with mirror slides out from either side of the table
- Integrated accessory rail for urinals, storage baskets etc.
- Shock-absorbent gas spring cylinder for easy height adjustment
- Two (2) Year Limited Warranty



### Specifications:

		table top Length		table top Width		Height Range		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
GF9000-LSTMO	Single Top, Laminate, Medium Oak, Graphite Grey U Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-LSTMC	Single Top, Laminate, Medium Cherry Graphite Grey U Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-LSTHM	Single Top, Laminate, Hardrock Maple, Graphite Grey U Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-LSTCL	Single Top, Laminate, Candlelight, Graphite Grey U Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-TSTMO	Single Top, Thermofoil, Medium Oak, Graphite Grey U-Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-TSTMC	Single Top, Thermofoil, Medium Cherry Graphite Grey U-Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-TSTHM	Single Top, Thermofoil, Hardrock Maple, Graphite Grey U-Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-TSTCL	Single Top, Thermofoil, Candlelight, Graphite Grey U-Base	32.5	82.55	19	48.26	28 - 44	71 - 111.8	34	15.42	100	45.36	1/ea
GF9000-TDTMO	Flip Top, Thermofoil, Medium Oak, Graphite Grey U-Base	32+13.5 (flip)=46	82.55+34.3=116.84	19	48.26	28 - 44	71 - 111.8	41	18.6	100	45.36	1/ea
GF9000-TDTMC	Flip Top, Thermofoil, Medium Cherry Graphite Grey U-Base	32+13.5 (flip)=46	82.55+34.3=116.84	19	48.26	28 - 44	71 - 111.8	41	18.6	100	45.36	1/ea
GF9000-TDTHM	Flip Top, Thermofoil, Hardrock Maple, Graphite Grey U-Base	32+13.5 (flip)=46	82.55+34.3=116.84	19	48.26	28 - 44	71 - 111.8	41	18.6	100	45.36	1/ea
GF9000-TDTCL	Flip Top, Thermofoil, Candlelight, Graphite Grey U-Base	32+13.5 (flip)=46	82.55+34.3=116.84	19	48.26	28 - 44	71 - 111.8	41	18.6	100	45.36	1/ea
GF9000-BSK	Storage Basket - 10X6x9 With Adapter											1/ea
GF9000-VNTY	Vanity Drawer With Mirror											1/ea
GF9000-UH	Urinal Hook With Adapter											1/ea
GF9000-CST	Non-Locking Caster											1/ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2019, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF-90000-9RevC22