

WHEELCHAIRS

TRANSIT TRANSPORT CHAIR

The EJT500 Transit Transport Chair is designed to safely and efficiently transport patients throughout the facility. Staff and caregivers can provide a comfortable seated transport experience while using a product designed to reduce bending and added strain during transport. The rigid frame is built to be a durable solution that also addresses the risk of facility theft.

Features:

- Easy to store with chair design that nests to save space when storing multiple chairs
- Weight capacity of 500 lb (226.8 kg) EVENLY DISTRIBUTED to accommodate use for various patients
- Foldable footrest for comfortable transport and fixed to avoid losing of parts
- Anti-tippers built into frame for added safety and stability
- Foldable armrest for easier entry and exit from chair
- Easy push handlebar with braking system
- Large under seat storage compartment
- Black Silvervein Hammertone Finish



Item	UOM	Overall Chair Width		Overall Chair Length		Overall Chair Height		Nested ¹ Length + each chair		Seat Width		Seat Depth		Seat Height from floor		Front Caster Diameter		Rear Wheel Diameter		Shipping Weight		Product Weight		Weight Capacity <small>EVENLY DISTRIBUTED</small>	
		inch [†]	cm	inch [†]	cm	inch [†]	cm	inch [†]	cm	inch [†]	cm	inch [†]	cm	inch [†]	cm	inch [†]	cm	inch [†]	cm	lb [†]	kg [†]	lb [†]	kg [†]	lb [†]	kg [†]
EJT5000	1/ea	26.77	68.00	41.929	106.50	40.354	102.50	12.00	30.48	17.717	45.00	18.11	46.00	21.654	55.00	5.00	12.70	8.00	20.32	70	31.75	65.2	29.57	500	226.80

[†] standard rounding method used

¹ approximate length footprint increases per chair width does not change



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website
www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF230008RevA23

Combined Footprint when nested for storage	2 Chairs		3 Chairs		4 Chairs		5 Chairs	
	inch [†]	cm	inch [†]	cm	inch [†]	cm	inch [†]	cm
Approximate Length width remains constant	54.00	137.16	66.00	167.64	78.00	198.12	90.00	228.60