



## PATIENT LIFTS AND SLINGS

## HYDRAULIC LIFT

Lumex floor lifts are designed to provide safer and more secure transfers of patients by caregivers. Lumex lifts are very sturdy and easy to use.

## **Features**

- Heavy gauge steel construction
- Six point spreader bar with 360° rotation
- Uses Lumex® six-point, four-point or two-point slings
- Foot pedal manually adjustable base
- Ergonomic placement of hydraulic pump handle facilitates ease of use
- Disassembles for easy storage
- Optional digital scale available
- Maximum weight capacity 400 lb - EVENLY DISTRIBUTED
- Designed to meet the requirments of HCPCS code:E0630

## Warranty

One (1) year limited warranty











|          |   | ov<br>Lei | erall<br>ngth | overall<br><b>Heiaht</b> |        | overall<br><b>Width</b> |       | Floor<br>Clearance |       | (meas | Spreader Bar Height (measured from sling attachment) minimum maximum |    |            | Base Width (outer measurement) minimum maximum |       |    |        | pro<br>We | oduct<br>eight | Weight<br>Capacity<br>EVENLY DISTRIBUTED |       |      |
|----------|---|-----------|---------------|--------------------------|--------|-------------------------|-------|--------------------|-------|-------|--|----|------------|--|-------|----|--------|-----------|----------------|--|-------|------|
|          |   | inch      |               | inch                     |        | _                       | cm    |                    |       | inch  |  |    | ich cm     |  |       | _  | cm     | _         |                | lb                                       |       |      |
| LF1030   | Blue Powder-coated                        | 48        | 121.92        | 55                       | 139.7  | 24.6                    | 62.48 | 5.5                | 13.97 | 27.4  | 69.6   | 80 | 0.3 203.96 | 24.6   | 62.48 | 43 | 109.22 | 90        | 40.8           | 400                                      | 181.4 | 1 ea |
| LF1031   | Black/Silver Vein Hammertone              | 48        | 121.92        | 55                       | 139.7  | 24.6                    | 62.48 | 5.5                | 13.97 | 27.4  | 69.6   | 80 | 0.3 203.96 | 24.6   | 62.48 | 43 | 109.22 | 90        | 40.8           | 400                                      | 181.4 | 1 ea |
| LF1031FP | Black/Silver Vein Hammertone w/foot pedal | 53.5      | 135.89        | 56.8                     | 144.27 | 26                      | 66.04 | 4.4                | 11.18 | 23.8  | 60.45  | 7  | 75 190.5   | 26   | 66.04 | 41 | 104.14 | 90        | 40.8           | 400                                      | 181.4 | 1 ea |



mGF Health Products, Inc. is an intertek ISO 13485:2016 / MDSAP Certified Company













