

LIFT ACCESSORIES

LIFT SCALE

Features

- Digital scale easily attaches to Lumex Patient Lifts
- Accurately weighs patients during daily lifts
- Combines the latest microprocessor and load cell technology to provide accurate and repeatable weight data
- Easy 4-button operation
- Weight is displayed for 5 minutes
- Auto power-off prolongs battery life
- Low battery indicator
- Measures pounds or kilograms
- Maximum weight capacity 600 lb
- Comes complete with attachment hardware

Warranty

Two (2) year limited warranty



		Weight Capacity EVENLY DISTRIBUTED		UOM
		lb	Kg	
DSC250	For use with LF1050 and LF1090	600	272.2	1 ea
DSC260	For use with LF1030/LF1031 and LF1040/LF1041	600	272.2	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF1200109RevC23