


LAMPS

DELUXE GOOSENECK EXAM LAMPS

Features

- Maximum of 60 Watt or 13 Watt CFL equivalent
- 115/120V AC. UL listed 3-wire plug with 8 ft (2.4 m) grounded cord
- Flared, doubled-wall, dual cool ventilated reflector
- Heavy gauge anodized aluminum gooseneck shade
- One piece phenolic socket with back-turn knob switch
- Height range: 34 ½" – 53 ½" (87.6 – 135.9 cm)



		minimum Height		maximum Height		UOM
		inch	cm	inch	cm	
1698	Chrome-plated Base, SOK-IT-GUARD Lock	34.5	87.6	53.5	135.9	1 ea
1698CM	Mobile Base/Clutch-Collar Lock	34.5	87.6	53.5	135.9	1 ea
1698M	Mobile Base, SOK-IT-GUARD Lock	34.5	87.6	53.5	135.9	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.