

DIAGNOSTIC INSTRUMENTS

GRAVES VAGINAL SPECULUM

Features

- Used for performing pelvic examinations, Pap smears, and examining the cervix
- It is a double blade duckbill speculum inserted with the blades closed
- Stainless steel



		overall Length		overall Width		UOM
		inch	cm	inch	cm	
2857	Small	3	7.6	0.75	1.9	1 ea
2858	Medium	4	10.2	1.25	3.2	1 ea
2859	Large	4.5	11.4	1.5	3.8	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.